



Item Number S1731
Item Sulfuric Acid, NF
CAS Number 7664-93-9
Molecular Formula H_2SO_4
Molecular Weight 98.08
Assay (by Wt.)..... 95.0-98.0%

MAXIMUM LIMITS

Residue on Ignition..... 0.005%
Chloride..... 0.005%
Arsenic..... 1 ppm
Heavy Metals..... 5 ppm
Reducing Substances..... To pass test
Residual solvents..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: East Coast Facility:
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843