



VWR ANALYTICAL

# Product Specification

**Catalog Number:** BDH1149-19L, BDH1149-4LG, 1, CABDH1149-19L, CABDH1149-4LG, CABDH1149-204L

**Catalog Description:** Tetrahydrofuran

**Grade:** ACS

**Date Current Specification Set:** January 15 2013

*All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.*

## Meets ACS Specifications for General Use

Characteristics	Specifications
Appearance:	Clear, colorless liquid
Water (K.F.):	Less than 0.050%
Residue:	Less than 0.030% <b>(0.0060/ 22.5mL)</b>
GC Purity:	Greater than 99.0%
Peroxides: solution)	Less than 0.015% as H <sub>2</sub> O <sub>2</sub> (No more than 3.9 mL of 0.1N sodium thiosulfate
Color (APHA):	Less than 20
Preservative (BHT): one cm cell wavelength, nm 278	Absorbance between 1.95 and 3.12 (250 - 400 ppm)

Signed on behalf of VWR:

Neal Fox  
Quality Control Supervisor