

# VWR glassware washer buying guide

GO TO  
VWR.COM

Please see  
VWR® Glassware  
Washer literature  
for complete  
specifications

## Standard washer or spindle rack-ready washer comparison

Standard Washer units with open racks:	Spindle Rack-Ready Washer units:
Ideal if your typical glassware includes beakers, test tubes, petri dishes, watch glasses and other wide-neck glassware.	Recommended for graduated cylinders, volumetric flasks, round bottom flasks and other narrow-neck glassware such as small Erlenmeyer flasks.
Includes 2 standard racks	Include 2 standard racks, can be customized with accessories and have DI Rinse function
Will not accept spindle racks	Designed to accept spindle racks which are sold separately
May have a DI rinse function	When spindle racks are used, additional plumbing connects the water feed directly to the spindle racks which inject it inside glassware
Employ spray arms which feed water underneath racks (similar to residential dish washers)	Can be used with: Full bottom spindle rack only (Cat. No. <b>10002-948</b> ) Full top and bottom spindle racks (Cat. No. <b>10002-950</b> ) Those with a combination of tall and short glassware can purchase bottom spindle rack and top-half spindle rack (Cat. No. <b>10002-946</b> )
Can be customized for various glassware with baskets and accessories	<b>Glassware Height:</b> 9 $\frac{3}{8}$ " clearance on the top spindle rack and 10 $\frac{1}{4}$ " clearance on the bottom rack; use bottom rack only for taller glassware – 20" maximum clearance

