




Claim Enhancements for Falcon® Serological Pipets

We are pleased to announce three claim enhancements to our Falcon serological pipet portfolio. In the spirit of continuously evolving to support the changing needs of our customers, we are implementing upgrades in the form of a human DNA-free claim, availability of extractables data, and increased sterility.



Human DNA-free

Falcon serological pipets will now feature a human DNA-free claim.

Pipets will be tested by DNA detection assay per ISO 18385 and found to be free of detectable human DNA.

Extractables data

Extractables analysis will be conducted in accordance with USP <665>. Data will be made available upon customer request.


Sterility Assurance Level (SAL) of 10⁻⁶

The serological pipets are sterilized and dosimetrically released per the requirements of ANSI/AAMI/ISO 11137.

The sterility of the upgraded products will be increased from SAL of 10⁻³ to 10⁻⁶.



Ordering Information



Falcon® Advantage® Packaging, Individually Wrapped

VWR Cat. No.	Corning Cat. No.	Size (mL)	Increments (mL)	Overall Capacity (mL)	Color Code	Qty/Pk	Qty/Cs
32314-012	356521	1	0.01	0.4	Yellow	100	1,000
32314-010	356507	2	0.01	0.8	Green	100	1,000
32314-008	356543	5	0.1	2.5	Blue	50	200
32314-006	356551	10	0.1	3.0	Red	50	200
32314-004	356525	25, Space Saver	0.25	7.0	Purple	50	200
32314-002	356535	25, Extended	0.5	11.0	Purple	50	200
32314-000	356550	50	1.0	10.0	Black	25	100
12777-018	357600	100	1.0	15.0	Black	5	50

Paper/Plastic Thermoform Packaging, Individually Wrapped

53300-283	357521	1	0.01	0.4	Yellow	100	1,000
53300-363	357507	2	0.01	0.8	Green	100	1,000
53300-421	357543	5	0.1	2.5	Blue	50	200
53300-523	357551	10	0.1	3.0	Red	50	200
53300-567	357525	25, Space Saver	0.25	7.0	Purple	50	200
53106-195	357535	25, Extended	0.5	11.0	Purple	50	200
53106-441	357550	50	1.0	10.0	Black	25	100

Bulk Packaged

53300-261	357506	1	0.01	0.4	Yellow	100	1,000
53300-365	357508	2	0.01	0.8	Green	100	1,000
53300-400	357529	5	0.1	2.5	Blue	50	200
53300-501	357530	10	0.1	3.0	Red	50	200
53392-198	357515	25, Space Saver	0.25	7.0	Purple	50	200

Literature Available

- ▶ Liquid Handling Ergonomics in Every Step
- ▶ Top 10 Ways to Reduce Contamination Risk
- ▶ Controlling Cell Culture Contamination

Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only.* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. *For a listing of US medical devices, regulatory classifications or specific information on claims, visit www.corning.com/resources.

Corning's products are not specifically designed and tested for diagnostic testing. Many Corning products, though not specific for diagnostic testing, can be used in the workflow and preparation of the test at the customers discretion. Customers may use these products to support their claims. We cannot make any claims or statements that our products are approved for diagnostic testing either directly or indirectly. The customer is responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of their intended application.

CORNING

For additional product or technical information, visit vwr.com/corning, call 1.800.932.5000, or contact your VWR representative.

 **avantor**[™]
delivered by **VWR**[™]

VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by country, federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. Offers valid in countries listed above, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2022 Avantor, Inc. All rights reserved.