



Certificate of Conformance, Analysis, and Sterility

This certification is provided as full assurance that the following product code and lot number were manufactured in accordance via a Quality System certified to ISO 9001 with prescribed procedures and specifications.

Product Code:	PSPS0076
Lot Number:	696327-0-0-1
Manufacturing Date:	04/06/2019
Expiration Date:	04/2021

Analysis

% v/v Isopropyl Alcohol:	70.46
Quantity of Particles Po:	5.6 millions/m ²
Quantity of Fibers Po:	0.25 thousands/m ²
Nonvolatile residue in DI:	0.00 g/m ²
Nonvolatile residue in IPA:	0.0205 g/m ²

Sterility

Gamma Processing Run ID:	2417862
Irradiation Dose Specified:	25 kGy to 50 kGy
Irradiation Dose Delivered:	29.1 kGy to 45.3 kGy
Validated Sterile Product:	Validation Sterile packets available upon request.
Sterility Assurance Level:	10 ⁻⁶ (SAL determined by test method AAMI/ISO 11137)

QC Supervisor
coc@contecinc.com

Thursday, April 18, 2019

NORTH AMERICA

Contec, Inc.
525 Locust Grove
Spartanburg, SC 29303 USA
tel: +1 864 503 8333
fax: +864 503 8444

EUROPE

Contec, Inc.
ZI du Prat - R.P. 3707
F-56037 VANNES Cedex France
tel: +33 (0)2 97 43 76 98
fax: +33 (0) 97 43 76 86

ASIA

Contec Cleanroom Technology (Suzhou) Co., Ltd.
17 Longyun Road
Suzhou Industrial Park, Suzhou 215024 China
tel: +86 512 6274 4050
fax: +86 512 6274 4051

Page 1 of 1

web: www.contecinc.com
e-mail: info@contecinc.com



Certificate of Processing

STERIGENICS 10811 Withers Cove Park Dr. Charlotte NC 28278
TEL 704 588-6877 FAX 704 588-8355 www.sterigenics.com

R55480102

04/15/2019 12:22:26 GMT
Page 1 of 1

Customer Name:	Contec- PO Box 530	Processing Facility:	Charlotte Gamma	Work Order #	2417862
P.O.#	99570			Sales Order #	2210263
	25-50 kGy	CNT1, Gamma Treatment	Irradiation Date/Time:	04/13/2019 07:37:00 GMT	

SO Line #	Qty	UOM	Customer Item Number	Customer Item Description	Customer Lot Number	Customer Load Number
101.000	54	CA	PSPS0076	Clean Room Product	696327-0-0-1	99570
102.000	54	CA	PSPS0076	Clean Room Product	696328-0-0-1	99570
103.000	54	CA	PSPS0076	Clean Room Product	696330-0-0-1	99570
104.000	108	CA	SB327030IR	32 oz IPA Spray Bottle 70/30	688950-0-0-1	99570
	270	CA	Total			

Quality Test Summary

Op#	Quality Test Description	Minimum Spec	Maximum Spec	Result	Pass/Fail	-----Signed By ----- User	Date /Time
450.00	Minimum Dose	25.0 kGy	50.0 kGy	29.1 KGY	Pass	MFREIBERG	04/13/2019 21:15:01 GMT
		Reason Code Test				MICHAEL FREIBERG	
450.00	Maximum Dose	25.0 kGy	50.0 kGy	45.3 KGY	Pass	BPOSTON	04/13/2019 19:36:37 GMT
		Reason Code Test				BRANTT POSTON	

Sterigenics certifies that the materials listed above (as described by the Manufacturer) received the indicated doses within the precision and accuracy of the dosimetry system employed.

Electronically Signed By: WENDY KOBLISH
Reason: Work Order Completions

Date: 04/15/2019 12:19:59 GMT