



Issuing Date 22-Jun-2010

Revision Date 10-May-2011

Revision Number 4

1. IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND THE COMPANY/UNDERTAKING

Product Name MES Free Acid Monohydrate

Product Code(s) E169

Recommended Use For Further Manufacturing Use Only
Not for Human or Animal Drug Use

Recommended Use Zwitterionic buffer for cell culture, enzyme assay, some electrophoretic applications

Synonyms MES (BUFFERING AGENT) * MORPHOLINOETHANESULFONIC ACID *
2-MORPHOLINOETHANESULFONIC ACID *

Company
AMRESCO, LLC
6681 Cochran Road
Solon, Ohio 44139

Emergency Telephone Number Chemtrec 1-800-424-9300

Company Phone Number 1-800-448-4442
E-mail Address info@amresco-inc.com

2. HAZARDS IDENTIFICATION

GHS - Classification

Classification Skin Irritant 2 Eye Irritant 2

Signal Word Warning



Hazard Statements H315 - Causes skin irritation H319 - Causes serious eye irritation H335 - May cause respiratory irritation

Precautionary Statements P280 - Wear protective gloves/protective clothing/eye protection/face protection
P302+ P352 - IF ON SKIN: Wash with plenty of soap and water P304 + P340 - IF
INHALED: Remove victim to fresh air and keep at rest in a position comfortable for
breathing P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for
several minutes. Remove contact lenses, if present and easy to do. Continue
rinsing

3. COMPOSITION/INFORMATION ON INGREDIENTS

Chemical Name	EC No.	REACH Reg. No.	CAS-No	Weight %	Classification
New Zealand assigned name pending	-	Not available	145224-94-8	95-100	-

4. FIRST AID MEASURES

Eye Contact Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician

Skin Contact Wash off immediately with soap and plenty of water removing all contaminated clothes and shoes

Inhalation Move to fresh air. If breathing becomes difficult, give oxygen

Ingestion Clean mouth with water and afterwards drink plenty of water

Notes to Physician Treat symptomatically

5. FIRE-FIGHTING MEASURES

Flammable Properties Not flammable

Flash Point Not determined

Suitable Extinguishing Media Dry chemical, CO2, water spray or regular foam

Hazardous Combustion Products Nitrogen oxides (NOx), Sulfur oxides, Carbon oxides.

Explosion Data

Sensitivity to Mechanical Impact Not sensitive

Sensitivity to Static Discharge Not sensitive

Protective Equipment and Precautions for Firefighters As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear

6. ACCIDENTAL RELEASE MEASURES

Personal Precautions Avoid contact with skin, eyes and clothing. Use personal protective equipment. Ensure adequate ventilation

Methods for Containment Prevent further leakage or spillage if safe to do so

Methods for Cleaning Up Avoid dust formation. Pick up and transfer to properly labeled containers. Ventilate area and wash spill site after material pickup is complete.

7. HANDLING AND STORAGE

Handling Handle in accordance with good industrial hygiene and safety practice.

Storage Keep containers tightly closed in a dry, cool and well-ventilated place

Incompatible Products Strong oxidizing agents.

8. EXPOSURE CONTROLS / PERSONAL PROTECTION

Chemical Name	ACGIH TLV	OSHA PEL	NIOSH IDLH
New Zealand assigned name pending 145224-94-8			

Engineering Measures Showers
Eyewash stations
Ventilation systems

Personal Protective Equipment

Eye/Face Protection Tightly fitting safety goggles

Skin and Body Protection Wear protective gloves/clothing

Respiratory Protection If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved respiratory protection should be worn. Positive-pressure supplied air respirators may be required for high airborne contaminant concentrations. Respiratory protection must be provided in accordance with current local regulations

Hygiene Measures Handle in accordance with good industrial hygiene and safety practice

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance	White to off-white
Odor	No information available
Odor Threshold	No information available
Decomposition Temperature °C	No information available
Freezing Point	No information available
Initial Boiling Point	No information available
Physical State	Crystalline Powder
pH	No data available
Flash Point	Not determined
Autoignition Temperature	No data available
Boiling Point/Range	No data available
Melting Point/Range	No data available
Flammability Limits in Air	Upper No data available Lower No data available
Explosive Properties	No information available
Oxidizing Properties	No information available
Evaporation Rate	No data available
MMHG @ 37.8 C	No data available
Vapor Density	No data available
Specific Gravity	No data available
Solubility	No information available
Partition Coefficient (n-octanol/water)	No data available
Viscosity	No information available

10. STABILITY AND REACTIVITY

Stability	Stable under recommended storage conditions
Incompatible Products	Strong oxidizing agents.
Conditions to Avoid	None known based on information supplied
Hazardous Decomposition Products	Carbon oxides. Nitrogen oxides (NOx). Sulfur oxides.
Hazardous Polymerization	Hazardous polymerization does not occur

11. TOXICOLOGICAL INFORMATION

Acute Toxicity

Product Information Irritating to eye, skin and respiratory system.

Chronic Toxicity

Chronic Toxicity No known effect based on information supplied

Target Organ Effects None known.

12. ECOLOGICAL INFORMATION

Ecotoxicity

The environmental impact of this product has not been fully investigated

13. DISPOSAL CONSIDERATIONS

Waste Disposal Method Dispose of material in accordance with all federal, state, and local regulations

Contaminated Packaging Dispose of in accordance with all federal, state and local regulations

14. TRANSPORT INFORMATION

IATA Not regulated

DOT Not regulated

15. REGULATORY INFORMATION

REACH Title VII Restrictions No information available.

Chemical Name	Dangerous Substances	Organic Solvents	Harmful Substances Whose Names Are to be Indicated on the Label	Pollution Release and Transfer Registry (Class II):	Pollution Release and Transfer Registry (Class I):	Poison and Deleterious Substances Control Law:
New Zealand assigned name pending	Not applicable	Not applicable	Not applicable	Not applicable	Not applicable	Not applicable

Component	ISHA - Harmful Substances Prohibited for Manufacturing, Importing, Transferring, or Supplying	ISHA - Harmful Substances Requiring Permission	Toxic Chemical Classification Listing (TCCL) - Toxic Chemicals	Toxic Release Inventory (TRI) - Group I	Toxic Release Inventory (TRI) - Group II
New Zealand assigned name pending 145224-94-8 (95-100)	Not applicable	Not applicable	Not applicable	Not applicable	Not applicable

International Inventories

TSCA	Does not Comply
EINECS/ELINCS	Does not Comply
ENCS	Does not Comply
IECSC	Does not Comply
KECL	Does not Comply
PICCS	Does not Comply
AICS	Does not Comply

U.S. Federal Regulations

SARA 313

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product contains a chemical or chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372

SARA 311/312 Hazard Categories

Acute Health Hazard	Yes
Chronic Health Hazard	No
Fire Hazard	No
Sudden Release of Pressure Hazard	No
Reactive Hazard	No

Clean Water Act

This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42)

Clean Air Act, Section 112 Hazardous Air Pollutants (HAPs) (see 40 CFR 61)

This product does not contain any substances regulated as hazardous air pollutants (HAPS) under Section 112 of the Clean Air Act Amendments of 1990

CERCLA

U.S. State Regulations

California Proposition 65

This product does not contain any Proposition 65 chemicals.

International Regulations

Mexico - Grade No information available.

Canada

WHMIS Hazard Class

D2B Toxic materials



16. OTHER INFORMATION

Issuing Date 22-Jun-2010

Revision Date 10-May-2011

Revision Note
No information available.

Recommended Restrictions No information available

Disclaimer

The information provided on this MSDS is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guide for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered as a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other material or in any process, unless specified in the text.

End of MSDS