



Product Specification

Catalog Number: BDH9242

Catalog Description: BDH L-ASCORBIC ACID ACS

Grade: ACS

Date Current Specification Set: NOV. 2012

All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.

Characteristics

Specifications

Heavy Metals (as Pb)	$\leq 0.002 \%$
Iron	$\leq 0.001 \%$
Purity	$\geq 99.0 \%$
Residue on Ignition	$\leq 0.1 \%$
Specific Rotation	+20.5 - +21.5 degrees

Signed on behalf of VWR:

Robin Crocker
Quality Control Manager