

**MATERIAL SAFETY DATA SHEET**

**Biotin Conjugated Antibody**

Rockland Immunochemicals, Inc.  
P.O. Box 326  
Gilbertsville, PA 19512

Emergency Telephone Number:  
Chemtrec (800) 424-9300

Information Telephone Number:  
(800) 656-7625

DATE PREPARED: 08-17-2006

REVIEWED:

**SECTION 1:-----CHEMICAL IDENTIFICATION-----**

NAME: BIOTIN CONJUGATED ANTIBODIES.

**SECTION 2:-----COMPOSITION/INFORMATION ON INGREDIENTS-----**

SUBSTANCE NAME CAS #

LIQUID ANTIBODY WITH .01% SODIUM AZIDE N/A

THE HAZARDS IDENTIFIED WITH THIS  
PRODUCT ARE THOSE ASSOCIATED WITH THE  
FOLLOWING COMPONENT(S):

SODIUM AZIDE 26628-22-8

SODIUM CHLORIDE 7647-14-5

**SECTION 3:-----HAZARDS IDENTIFICATION-----**

EMERGENCY OVERVIEW

SODIUM AZIDE MAY REACT WITH LEAD AND COPPER PLUMBING TO FORM  
HIGHLY EXPLOSIVE METAL AZIDES.

HMS RATING

HEALTH: 0

FLAMMABILITY: 0

REACTIVITY: 1

NFPA RATING

HEALTH: 0

FLAMMABILITY: 0

REACTIVITY: 1

FOR ADDITIONAL INFORMATION ON TOXICITY, PLEASE REFER TO SECTION 11.

**SECTION 4:-----FIRST-AID MEASURES-----**

ORAL EXPOSURE

IF SWALLOWED, WASH OUT MOUTH WITH WATER PROVIDED PERSON IS  
CONSCIOUS. CALL A PHYSICIAN IMMEDIATELY.

INHALATION EXPOSURE

IF INHALED, REMOVE TO FRESH AIR. IF NOT BREATHING GIVE ARTIFICIAL  
RESPIRATION. IF BREATHING IS DIFFICULT, GIVE OXYGEN.

DERMAL EXPOSURE

IN CASE OF SKIN CONTACT, FLUSH WITH COPIOUS AMOUNTS OF WATER FOR AT  
LEAST 15 MINUTES. REMOVE CONTAMINATED CLOTHING AND SHOES. CALL A  
PHYSICIAN.

EYE EXPOSURE

In case of contact with eyes, flush with copious amounts of water for at least 15 minutes. Assure adequate  
flushing by separating the eyelids with fingers. Call a physician.

**SECTION 5:-----FIRE FIGHTING MEASURES-----**

**EXPLOSION HAZARDS**

AZIDE REACTS WITH MANY HEAVY METALS SUCH AS LEAD, COPPER, MERCURY, SILVER, GOLD TO FORM EXPLOSIVE COMPOUNDS. COPPER AND LEAD AZIDES ARE MORE SENSITIVE THAN NITROGLYCERINE. AZIDE REACTS WITH METAL HALIDES TO GIVE A RANGE OF METAL AZIDE HALIDES, MANY OF WHICH ARE EXPLOSIVE. INCOMPATIBLE WITH CHROMYL CHLORIDE, HYDRAZINE, BROMINE, CARBON DISULFIDE, DIMETHYL SULFATE, DIBROMOMALONITRILE.

**AUTOIGNITION TEMP**

N/A

**FLAMMABILITY**

N/A

**EXTINGUISHING MEDIA**

SUITABLE: WATER SPRAY, CARBON DIOXIDE, DRY CHEMICAL POWDER, OR APPROPRIATE FOAM.

PROTECTIVE EQUIPMENT: WEAR SELF-CONTAINED BREATHING APPARATUS AND PROTECTIVE CLOTHING TO PREVENT CONTACT WITH SKIN AND EYES.

SPECIFIC HAZARD(S): EMITS TOXIC FUMES UNDER FIRE CONDITIONS.

**SECTION 6:-----ACCIDENTAL RELEASE MEASURES-----**

**PROCEDURE(S) OF PERSONAL PRECAUTION(S)**

WEAR RESPIRATOR, CHEMICAL SAFETY GOGGLES, RUBBER BOOTS, AND HEAVY RUBBER GLOVES.

**METHODS FOR CLEANING UP:**

SPILLED MATERIAL SHOULD BE CAREFULLY WIPED UP OR MOISTENED WITH WATER AND REMOVED. VENTILATE AREA AND WASH SPILL SITE AFTER MATERIAL PICKUP IS COMPLETE.

**SECTION 7:-----HANDLING AND STORAGE-----**

**HANDLING**

USER EXPOSURE: AVOID INHALATION. AVOID CONTACT WITH EYES, SKIN, AND CLOTHING. AVOID PROLONGED OR REPEATED EXPOSURE.

**SECTION 8:-----EXPOSURE CONTROLS/PERSONAL PROTECTION-----**

**PERSONAL PROTECTIVE EQUIPMENT**

RESPIRATORY: USE RESPIRATORS AND COMPONENTS TESTED AND APPROVED UNDER APPROPRIATE GOVERNMENT STANDARDS SUCH AS NIOSH (US) OR CEN (EU). WHERE RISK ASSESSMENT SHOWS AIR-PURIFYING RESPIRATORS ARE APPROPRIATE USE A FULL-FACE PARTICLE RESPIRATOR TYPE N100 (US) OR TYPE P3 (EN 143) RESPIRATOR CARTRIDGES AS A BACKUP TO ENGINEERING CONTROLS. IF THE RESPIRATOR IS THE SOLE MEANS OF PROTECTION, USE A FULL-FACE SUPPLIED AIR RESPIRATOR.

HAND: COMPATIBLE CHEMICAL-RESISTANT GLOVES.

EYE: CHEMICAL SAFETY GOGGLES.

**GENERAL HYGIENE MEASURES**

WASH THOROUGHLY AFTER HANDLING. WASH CONTAMINATED CLOTHING BEFORE REUSE.

**SECTION 9:-----PHYSICAL AND CHEMICAL PROPERTIES-----**

**APPEARANCE AND ODOR**

LIQUID

**SECTION 10:-----STABILITY AND REACTIVITY -----**

**STABLE**

MATERIALS TO AVOID: DIMETHYL SULFATE IS INCOMPATIBLE WITH SODIUM AZIDE, ACID CHLORIDES, HALOGENATED SOLVENTS AVOID CONTACT WITH METALS. AVOID CONTACT WITH ACID. SODIUM AZIDE MAY REACT WITH LEAD AND COPPER PLUMBING TO FORM HIGHLY EXPLOSIVE METAL AZIDES.

SODIUM CHLORIDE: AVOID STRONG OXIDIZING AGENTS AND ACIDS

HAZARDOUS DECOMPOSITION PRODUCTS: NATURE OF DECOMPOSITION

PRODUCTS IS UNKNOWN

HAZARDOUS POLYMERIZATION: WILL NOT OCCUR

**SECTION 11:-----TOXICOLOGICAL INFORMATION-----**

ROUTE OF EXPOSURE

SKIN CONTACT: MAY CAUSE SKIN IRRITATION.

SKIN ABSORPTION: MAY BE HARMFUL IF ABSORBED THROUGH THE SKIN.

EYE CONTACT: MAY CAUSE EYE IRRITATION.

INHALATION: MAY BE HARMFUL IF INHALED. MATERIAL MAY BE IRRITATING TO MUCOUS MEMBRANES AND UPPER RESPIRATORY TRACT.

INGESTION: MAY BE HARMFUL IF SWALLOWED.

SIGNS AND SYMPTOMS OF EXPOSURE:

MANY AZIDES CAUSE A FALL IN BLOOD PRESSURE AND SOME INHIBIT ENZYME ACTION. LABORATORY EXPERIMENTS IN ANIMALS HAVE SHOWN SODIUM AZIDE TO PRODUCE A PROFOUND HYPOTENSIVE EFFECT, DEMYELINATION OF MYELINATED NERVE FIBERS IN THE CENTRAL NERVOUS SYSTEM, TESTICULAR DAMAGE, BLINDNESS, ATTACKS OF RIGIDITY, AND HEPATIC AND CEREBRAL EFFECTS.

TO THE BEST OF OUR KNOWLEDGE, THE CHEMICAL, PHYSICAL, AND TOXICOLOGICAL PROPERTIES HAVE NOT BEEN THOROUGHLY INVESTIGATED.

**SECTION 12:-----ECOLOGICAL INFORMATION-----**

DATA NOT YET AVAILABLE

**SECTION 13:-----DISPOSAL CONSIDERATIONS-----**

APPROPRIATE METHOD OF DISPOSAL OF SUBSTANCE OR PREPARATION

CONTACT A LICENSED PROFESSIONAL WASTE DISPOSAL SERVICE TO DISPOSE OF THIS MATERIAL. DISSOLVE OR MIX THE MATERIAL WITH A COMBUSTIBLE SOLVENT AND BURN IN A CHEMICAL INCINERATOR EQUIPPED WITH AN AFTERBURNER AND SCRUBBER. OBSERVE ALL FEDERAL, STATE, AND LOCAL ENVIRONMENTAL REGULATIONS.

**SECTION 14:-----TRANSPORT INFORMATION-----**

DOT

PROPER SHIPPING NAME: NONE

NON-HAZARDOUS FOR TRANSPORT: THIS SUBSTANCE IS CONSIDERED TO BE NON-HAZARDOUS FOR TRANSPORT.

IATA

NON-HAZARDOUS FOR AIR TRANSPORT: NON-HAZARDOUS FOR AIR TRANSPORT.

**SECTION 15:-----REGULATORY INFORMATION-----**

NOT AVAILABLE.

**SECTION 16:-----OTHER INFORMATION-----**

THE ABOVE INFORMATION IS BELIEVED TO BE CORRECT BUT DOES NOT PURPORT TO BE ALL INCLUSIVE AND SHALL BE USED ONLY AS A GUIDE. MANUFACTURER SHALL NOT BE HELD LIABLE FOR ANY DAMAGE RESULTING FROM HANDLING FROM CONTACT WITH THE ABOVE PRODUCT.

**CHECK OUT OUR:**

[PHOSPHO SPECIFIC ANTIBODIES](#)

[MOLECULAR BIOLOGY BUFFERS](#)

[CUSTOM ANTIBODY SERVICES](#)