



Material	82120-300
Material description	Sodium cyanide
Grade	Reagent
CAS Number	143-33-9
Molecular formula	CNNa
Molecular mass	49,01

Characteristics	Specifications
-----------------	----------------

Assay	Min. 95,0 %
Total Fe (Iron)	Max. 0,005 %
Cl (Chloride)	Max. 0,15 %
PO4 (Phosphate)	Max. 0,05 %
S (Sulphide)	Max. 0,005 %
SCN (Thiocyanate)	Max. 0,02 %
SO4 (Sulphate)	Max. 0,05 %
Pb (Lead)	Max. 5 ppm

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Michel Boisvert, Head of the Laboratory - Lachine
VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3,
QC, Canada