

Ovation® / M

Mechanical Volume-Set Pipettes




Performance, ergonomics and value in one elegant design.

Actual Size



Simple conventional volume setting



Easy tip acquisition
EVERY TIME




Ergonomically-correct posture that ensures pipetting comfort

Innovation you can feel.



Fits your hand like a glove.

Ovation's contoured shape and adjustable hook means a custom fit for every user in the lab.



Easy tip acquisition & discard.

An audible "click" indicates tips are securely seated and ready for use. And discarding tips takes just the gentle push of a button.



Light thumb action.

A very low-force plunger allows the user to manually control the speed of aspirating and dispensing with minimal effort.



Fully Autoclavable

No disassembly needed.




**6 autoclavable models
to cover 0.2–1000µL**



Conveniently open PCR and microtubes.

Ovation has a cap opener built right into its base.



Simple in-lab calibration.

No need to send the pipette out.

Order Information

VWR Catalog No.	Color	Volume Range	Dispensing Increments	Accuracy*	Precision*
10018-836	Red	0.2 – 2µL	0.002µL	±1.5% at 2µL ±2.7% at 1µL ±12% at 0.2µL	0.9% at 2µL 1.8% at 1µL 7.6% at 0.2µL
10018-838	Orange	1 – 10µL	0.01µL	±1% at 10µL ±1.5% at 65µL ±2.5% at 2µL	0.4% at 10µL 0.6% at 5µL 1.2% at 1µL
10018-840	Yellow	2 – 20µL	0.02µL	±1% at 20µL ±1.5% at 10µL ±7% at 2µL	0.3% at 20µL 0.5% at 10µL 2% at 2µL
10018-842	Green	10 – 100µL	.1µL	±0.8% at 100µL ±0.8% at 50µL ±3.5% at 10µL	0.15% at 100µL 0.24% at 50µL 1% at 10µL
10018-844	Blue	20 – 200µL	.2µL	±0.8% at 200µL ±0.8% at 100µL ±2.5% at 20µL	0.15% at 200µL 0.25% at 100µL 1% at 20µL
10018-846	Purple	100 – 1000µL	1.0µL	±0.8% at 1000µL ±0.8% at 500µL ±3% at 100µL	0.15% at 1000µL 0.2% at 500µL 0.6% at 100µL

* Specifications subject to change



VistaLab tips – A fitting companion.

Whether the requirements of your application require tips that are sterile, certified RNase/DNase free, non-pyrogenic, Trace Metal certified or filtered, there's a genuine VistaLab tip to meet your needs.

Ovation M Starter Sets

3 and 4 pipette starter sets are available at special pricing.

Order **10018-848** for the 1-10µL, 10-100µL and 100-1000µL Set

Order **10018-850** for the 0.2-2µL, 2-20µL and 20-200µL and 100-1000µL Set

VWR  **1.800.932.5000**
We Enable Science **vwr.com**

Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2016 VWR International, LLC. All rights reserved.