



Material	58934-940
Material description	Methyl ethyl ketone
Grade	ACS
CAS Number	78-93-3
Molecular formula	C ₂ H ₅ COCH ₃
Molecular mass	72.11

Characteristics

Specifications

Assay (GC, corrected for water)	Min. 99.0 %
Titration acid	Max. 0.0005 meq/g
Color (Pt-Co)	Max. 15 APHA
Residue after evaporation	Max. 0.0025 %
Water	Max. 0.20 %

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Michel Boisvert, Head of the Laboratory - Lachine
VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3,
QC, Canada

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.