



## VWR® BASIC PROTECTION SPP SHOE COVERS

As per the manufacturer, the below product meets the following criteria:

<b>VWR Catalog Number:</b>	414004-653
<b>Lot Number:</b>	31920193
<b>Description:</b>	VWR Shoecvr Antskd Blu Xlcs300.
<b>Date of Manufacture:</b>	November 22, 2019
<b>Country of Origin:</b>	CHINA

### Quality System Compliance

Products are manufactured under the ISO13485:2003 standard. Products are Inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP. Inspection records are reviewed and signed off by qualified personnel for product release

### QC Testing

Inspected according to AQL level II 4.0 sampling plan per MIL-STD-105E

<b>Break Tensile Properties (N):</b>	GB/T3923.1-1997 standard	MD = 48.55	CD = 36.83
<b>Elongation Rate (%):</b>	GB/T3923.1-1997 standard	MD = 78.93	CD = 79.32

### Product Specifications

This product is a Class I category device as defined by the FDA in 21CFR Parts 862-892

<b>Material:</b>	35gsm SPP-Blue (Spunbonded polypropylene fabric)
<b>BSE/STE Statement:</b>	The product contains no ingredients of animal origin and no material derived from or exposed to animals affected by, or under quarantine, for transmitting animal spongiform encephalopathy/bovine spongiform encephalopathy.
<b>Non-Pyrogenic:</b>	N/A
<b>DNase &amp; RNase Free:</b>	N/A
<b>Cytotoxicity:</b>	N/A
<b>Latex Free Statement:</b>	This product is latex free.

Signed:



Jamie Ethier  
VP Global Quality  
VWR, Part of Avantor

Date: October 19, 2020