

**Eldon James****Certificate of Conformance/Traceability**

**Customer Name:** VWR International\*  
**P.O. Number:** 4508935881

**Pack Slip Number:** 164000  
**Ship Date:** 9/12/18

**Part Number:** C8-6NK-QC

**Customer Part #:** 89238-118

**Lot Number:** 1232792

**Customer Rev :**

**Quantity:** 200

**Manufacture Date:** 6/12/2018 12:00:00AM

**Raw Material:** Natural Kynar® - Polyvinylidene Fluoride (PVDF)

**Material Grade:** Kynar® 1000HD

**Manufacturer:** Arkema AKA aka Elf Atochem aka Atofina

These products are manufactured in the United States of America by the Eldon James Corp. Eldon James, Corp. fittings with a suffix of NK are molded from Virgin ingredients supplied by the manufacturer listed above. Eldon James, Corp. certifies that all fittings meet the specifications as stated in current catalogs and current Quality Control drawings. Eldon James, Corp. will not take responsibility for errors or changes customers make to in-house drawings of Eldon James, Corp. products. Customers are urged to request samples if they are uncertain of what they are purchasing. Ultimately, customers must make their own determination that their use of our product(s) is safe, lawful (except as provided in the above statements), and technically suitable for their intended application(s).

**Kynar 1000HD Raw Material Regulatory Approvals:** NSF-61; NSF 14; FDA 21 CFR 177.2510 (repeated contact with food) FDA 21 CFR 177.1520 (1% concentration single use adjuvant for Polyolefins). This material is not known to contain latex, phthalates, melamine, GMOs or VOCs (volatile organic compounds).


**US Pharmacopoeia (USP) Kynar®** This material has been certified to the **USP Class VI Plastic (USP<88>**, Biological Reactivity Tests, In Vivo) requirements. Certificate is available upon request.


**Gamma** This material is gamma irradiation stable.

**Tallow/BSE/TSE** Based on information supplied to us by the resin manufacturer, Kynar 1000HD is Animal Derivative Free (ADF) "We do not use tallow or any other animal derived materials in the manufacture or formulation of this product".

**RoHs/REACH** This product is RoHs compliant. Based on information provided by Arkema, Kynar 1000HD contains no SVHCs (Substances of Very High Concern) and meets REACH requirements. No additional additives are employed in the manufacture of Eldon James Natural Kynar (NK) fittings.

For Additional information about the resin material used in this product, please request a PDF copy of the manufacturer's **Regulatory Information/Product Stewardship Data Sheet** from the Eldon James Corp.

  
 \_\_\_\_\_  
 Eldon James Corp.

  
 \_\_\_\_\_  
 Date: