




AntiStatic Kit Discharge Your Samples




Pass your sample through the U-ionisator
The AntiStatic Kit completely removes all electrostatic charge from your samples. Scattering of powdery weighing substances is eliminated, ensuring highly reproducible results.



AntiStatic Kit

Equip your XP Analytical, Micro and Ultra Micro Balances with an AntiStatic Kit. The U-Electrode discharges your sample automatically before weighing – absolutely safe to the touch.


AntiStatic Kits



Part **Universal AntiStatic Kit**
complete (incl. U-Electrode large and Universal Power Supply)

Balances all balances

Order No. 11107767



Integrable AntiStatic Kit
complete (incl. Point Electrode and XP Power Supply)

XP Analytical Balances and XP56/26 (incl. DR) Micro Balances

Order No. 11107761

Second Point Electrode
optional for XP integrable AntiStatic Kit

Order No. 11107762

Single AntiStatic Components



Part **U-Electrode large**

Recommended for discharging all kind of samples and tare containers

Balances all balances

Order No. 11107764



U-Electrode small

Recommended for discharging small samples

Balances all balances

Order No. 11140161



Point Electrode

Recommended for discharging small samples on spatulas

Balances all balances

Order No. 11107765

XP Power Supply

Automated switch on / off, synchronized with door opening

Balances XP Analytical, Micro and Ultra Micro Balances

Order No. 11107763

Universal Power Supply

Manual switch on / off

Balances all balances

Order No. 11107766

Mettler-Toledo AG

Laboratory & Weighing Technologies
CH-8606 Greifensee, Switzerland
Phone +41 44 944 22 11
Fax +41 44 944 30 60

Subject to technical changes
© 01/2010 Mettler-Toledo AG
Printed in Switzerland 1 1796356
Global MarCom Switzerland

www.mt.com/lab-AntiStatic

For more information