


mdi

Membrane Technologies

Stericheck

Closed Sterility Test Systems



- Minimizes risk of false positives/negatives
- Meets global regulatory requirements

Pre-sterilized, Ready to use

- Gamma irradiated
- EO Sterilized



Validated for

- Sterility
- Microbial Recovery
- Microbial Retention
- Equal Sample Distribution



Types Available

Stericheck SVP1
Small Volume Parenterals
Liquid Vials

Stericheck SVP2
Ampules & Collapsible bags

Stericheck SVP3
Small Volume Parenterals
Vials with soluble powder

Stericheck SVP4
Small Volume Parenterals
Mini Vials and Cartridges

Stericheck PC
Plastic containers

Stericheck LVP
Large Volume Parenterals
in glass bottles


Stericheck PFS
Pre-filled Syringes

Stericheck MPFS
For Medical Devices with
Female Luer Connection

Stericheck ML
For Medical Devices with
Male Luer Connection

SVP Dilutor
Antibiotics and difficult to
dissolve powders

Unique Performance Advantage



Non Coring Needle Design

- To prevent blockage by rubber particles from vial bungs

Pre-installed Color Coded Clamps

- For easy identification and clamping at both ends of tube




Raised Transparent Vent

- To minimize chances of fluid logging and subsequent obstruction of air flow
- For easy visibility in case of accidental fluid logging

Special 'L' Shaped Flow Directors

- To minimize frothing, specially with viscous fluids



Tough Membrane

- Withstand inadvertent back pressure
- Available with three different polymeric tough membranes

Archivable Chemical Indicator

- In each Stericheck pack as evidence of EO gas sterilization



Easy to Peel off, Impervious Blister Pack Cover

- For radiation sterilized Stericheck devices to prevent ingress of H₂O₂ during isolator disinfection

Customized Needle Connections

- For specific customer needs such as mini vials and cartridges in biopharmaceuticals



Ordering Information

Product Type	Membrane (0.45 µm)	EO Sterile		Gamma Sterile	
		MDI Catalogue no.	VWR Cat. no.	MDI Catalogue no.	VWR Cat. no.
Stericheck SVP1	PVDF	SV1V0902XXXX202	75792-168	SV1V0902XXXX302	75792-172
	CN	SV1C0902XXXX202	75792-170	SV1C0902XXXX302	75792-174
Stericheck SVP2	PVDF	SV2V0902XXXX202	75792-176	SV2V0902XXXX302	75792-180
	CN	SV2C0902XXXX202	75792-178	SV2C0902XXXX302	75792-182
Stericheck SVP3	PVDF	SV3V0902XXXX202	75792-184	SV3V0902XXXX302	75792-044
	CN	SV3C0902XXXX202	75792-186	SV3C0902XXXX302	75792-046
Stericheck SVP4	PVDF	SV4V0902XXXX202	75792-048	SV4V0902XXXX302	75792-052
	CN	SV4C0902XXXX202	75792-050	SV4C0902XXXX302	75792-054
Stericheck LVP	PVDF	SLVV0902XXXX202	75792-056	SLVV0902XXXX302	75792-090
	CN	SLVC0902XXXX202	75792-058	SLVC0902XXXX302	75792-092
Stericheck PC	PVDF	SPCV0902XXXX202	75792-094	SPCV0902XXXX302	75792-060
	CN	SPCC0902XXXX202	75792-096	SPCC0902XXXX302	75792-062
Stericheck PFS	PVDF	SPFV0902XXXX202	75792-064	SPFV0902XXXX302	75792-098
	CN	SPFC0902XXXX202	75792-066	SPFC0902XXXX302	75792-100
Stericheck Dilutor	-	SVPD0902XXXX209	75792-102	SVPD0902XXXX309	75792-104
Stericheck MPFS	PVDF	SMPV090201XX202	75792-106	SMPV090201XX302	75792-110
	CN	SMPC090201XX202	75792-108	SMPC090201XX302	75792-112
Stericheck ML	PVDF	SMLV0902XXXX202	75798-778	SMLV0902XXXX302	75796-902
	CN	SMLC0902XXXX202	75798-780	SMLC0902XXXX302	75796-904

VWR  **1.800.932.5000**
We Enable Science **vwr.com**

Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2017 VWR International, LLC. All rights reserved.