



Item Number S1525

Item Sorbitol Solution, 70%, FCC

CAS Number

Molecular Formula

Molecular Weight

Assay..... Min. 64.0%

pH of a 14% (w/w) Solution..... 5.0-7.5

Water.....28.5-31.5%

MAXIMUM LIMITS

Lead..... 1 mg/kg

Nickel..... 1 mg/kg

Reducing Sugars.....0.3%

Residue on Ignition..... 0.1%

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: East Coast Facility:

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843