

Edetate Disodium
U.S.P.
(edta, disodium salt dihydrate crystal)



Material No.: 8994-01
Revision No.: 0

Product Specification

Meets U.S.P Requirements, GMP Manufactured Product

Test	Specification
USP – Identification A	Passes Test
USP – Identification B	Passes Test
USP – Identification C	Passes Test
USP – pH (1 in 20)	4.0 – 6.0
USP – Loss on Drying	8.7 – 11.4 %
USP – Calcium (Ca)	Passes Test
USP – Nitritotriacetic Acid	≤ 0.1 %
USP – Assay (C ₁₀ H ₁₄ N ₂ Na ₂ O ₈) (dried basis)	99.0 – 101.0 %

GMP Manufactured Product
Bulk Pharmaceutical Chemical

CAUTION: For Manufacturing, processing or repackaging

No Class 1,2,3 or other solvents are used or produced in the manufacturing or purification of the product.

Elemental Impurities (USP 232, EP 5.20) – Information on elemental impurities for this product is available on the associated Product Regulatory Data Sheet and elemental impurity profile report.

Packaging Site: Paris Mfg Ctr & DC


Jamie Ethier
Vice President Global Quality

For questions on this Product Specification please contact Technical Services at 855.282.6867 or +1.610.386.1700

Avantor Performance Materials, LLC

100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone 610.386.1700

Page 1 of 1