



Item Number A1178
Item Ammonium Dichromate, Crystal, Reagent
CAS Number 7789-09-5
Molecular Formula $(\text{NH}_4)_2\text{Cr}_2\text{O}_7$
Molecular Weight 252.10

Assay..... Min. 99.5%
pH of a 5% Solution @ 25°C..... 3-5

MAXIMUM LIMITS

Alkali (as SO_4).....0.1%
Calcium (Ca)..... 0.005%
Chloride (Cl)..... 0.005%
Insoluble Matter..... 0.005%
Iron (Fe).....0.005%
Sulfate (SO_4).....0.01%

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: 14422 S. San Pedro St.
Gardena, CA 90248
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843

East Coast Facility: 755 Jersey Ave.
New Brunswick, NJ 08901