



VWR® MAXIMUM PROTECTION FLUID-IMPERVIOUS BUTTERFLY-STYLE SHOE COVERS

As per the manufacturer, the below product meets the following criteria:

VWR Catalog Number:	414004-643
Lot Number:	03192093
Description:	VWR SHOECVR FLUID BLUE 2INX 75PR
Date of Manufacture:	December 10, 2019
Country of Origin:	CHINA

Quality System Compliance

Products are manufactured under the ISO 9001:2015 & ISO13485:2003 standard. Products are Inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP. Inspection records are reviewed and signed off by qualified personnel for product release.

QC Testing

Inspected according to AQL level II 4.0 sampling plan per MIL-STD-105E

-Break Tensile properties(N):	ASTM D5034standard,	MD= 99.62	CD= 105.86
-Elongation rate (%):	ASTM D5034standard,	MD= 72.78	CD= 88.41
-Trapezoid tearing strength(N):	FZ/T 01085-1999 standard,	MD= 639.00	CD= 631.00

Product Specifications

Material:	SPP+PE 60g/m ² Blue (Spunbonded Polypropylene fabric+polyethylene fabric)
BSE/TSE Statement:	The above product contains no ingredients of animal origin and no material derived from, or exposed to animals affected by or under quarantine for Transmitting, Animal Spongiform Encephalopathy (TSE)/Bovine Spongiform Encephalopathy(BSE)
Non-Pyrogenic:	NA
DNase & RNase Free:	NA
Cytotoxicity:	NA
Latex Free Statement:	The above product is latex free.

Signed:



Jamie Ethier
VP Global Quality
VWR, Part of Avantor

Date: January 04, 2020