



BIOFIT SEATING FOR
HEALTHCARE &
LABORATORY
ENVIRONMENTS
**FIGHTING
INFECTION –
THE BIOFIT
ADVANTAGE**



PROTECTING THE CAREGIVER AND BEYOND

Now more than ever, easy, effective cleaning and disinfection of healthcare and laboratory environments, equipment and furnishings are critical. That's why BioFit developed a comprehensive collection of ergonomic seating specially designed with durable materials that stand up to more frequent and longer exposure to cleaners and disinfectants.

VWR carries numerous seating models from the BioFit Quick Ship program that are available for next-day or 2-day shipping, including ISO 5 and ISO 6 cleanroom configurations, as well as the other designs featured in this catalog. Of note, the BioFit Rexford stool included in the fighting infection product lineup is ideal for emerging healthcare needs and available for a limited period at a reduced lead time of 5 business days.





Tight-fitting, seamless plastic backs ensure a smooth, easy-to-clean surface.

Non-porous, coated, chemical-resistant vinyl used on the seats and backs can be cleaned and disinfected with EPA certified cleaning and disinfecting products. Durable coating and BioFit's 13-year warranty means they'll stand up to frequent cleaning and extended disinfectant exposure.

Smooth plastic controls wipe down easily.

Seamless welded footrings with chrome or black powdercoat finishes are easy to clean and disinfect because of their smooth surface.

Footring supports feature the same T-shaped profile as the base, with no recesses for superior cleaning and disinfecting.

5-Star Cast Aluminum or Resin Bases offer great stability and feature a T-shaped profile for easy cleaning and disinfecting. Models with steel bases and affixed footrings feature seamless welded tubing and chrome or black powder coat finishes.

Hooded resistance casters feature smooth tops and sides for easy cleaning and disinfecting.

BIOFIT QUICK SHIP SEATING – WHEN YOU NEED CHAIRS FAST*

Need chairs in a hurry? BioFit offers next-day shipping (NDS) and 2-day shipping (2DS) in limited quantities on a variety of popular configurations for lab and healthcare applications – indicated by an orange dot ● next to eligible NDS models – and a green dot ● next to eligible 2DS models.

ELITE SERIES SEATING

Available with non-porous PVC chemical-resistant vinyl upholstery which can be cleaned and disinfected.

- This versatile series meets the needs of nearly any healthcare facility workspace or lab discipline, including clean rooms rated ISO 5 and ISO 6.
- Roomy 18.5” wide x 17” deep x 3” thick seat with four-way contouring and internal seat board bumperguard and lumbar-supporting 15.5” wide x 14.5” high backrest provide all-day comfort.
- Full ergonomic functionality and seat height ranges to accommodate many tasks.



VWR Catalog #	BioFit Part #	Seat Height	Control Style	Foot Type	Upholstery	Footring	Base Style	Metal Finish	Price
● 76404-332	EEAMHFFXXC540	19-26"	Free Float Articulating	Glides	Black Vinyl	20" Adjustable	High Profile Aluminum	Chrome	\$1,083.40
● 76404-334	EESHHTXXC540	25-32"	Task	Glides	Black Vinyl	21" Affixed	Tubular Steel	Chrome	\$1,023.26

ELITE SERIES ERGONOMIC CHAIR WITH HEIGHT-ADJUSTABLE ARMS

VWR Catalog #	BioFit Part #	Seat Height	Control Style	Foot Type	Upholstery	Footring	Base Style	Metal Finish	Price
● 76404-354	EERMRFPPXC540	19-26"	Free Float Articulating	Casters	Black Vinyl	NA	Reinforced Composite	Chrome	\$1,080.97

ELITE SERIES ISO 5 CLEAN ROOM CHAIR

VWR Catalog #	BioFit Part #	Seat Height	Control Style	Foot Type	Upholstery	Footring	Base Style	Metal Finish	Price
● 76404-360	EESMHTXX5C540	21-28"	Task	Glides	Black Vinyl	21" Affixed	Tubular Steel	Chrome	\$1,053.33

ELITE SERIES ISO 6 CLEAN ROOM CHAIRS

VWR Catalog #	BioFit Part #	Seat Height	Control Style	Foot Type	Upholstery	Footring	Base Style	Metal Finish	Price
● 76404-356	EESHHTXX6C540	25-32"	Task	Glides	Black Vinyl	21" Affixed	Tubular Steel	Chrome	\$912.64
● 76404-362	EESMRTXX6C540	21-28"	Task	Casters	Black Vinyl	21" Affixed	Tubular Steel	Chrome	\$900.62

*Based on model, Quick Ship products ship next business day (NDS) or in 2 business days (2DS) in limited quantities per order (9 units per NDS order; 15 models per 2DS order).

UniqueU SERIES URETHANE SEATING

Features non-porous urethane seat and back which can be cleaned and disinfected; resistant to water, many chemicals, and punctures.

- These durable chairs are engineered for whatever the workday can dish out.
- Comfy 18.5" wide x 17" deep x 1.5" thick seat, lumber-supporting 16.5" wide x 12.5" high backrest and ergonomic functionality combine to deliver comfort and performance.
- Available with pneumatic seat height adjustment ranges to accommodate many tasks.



VWR Catalog #	BioFit Part #	Seat Height	Control Style	Foot Type	Upholstery	Footring	Base Style	Metal Finish	Price
• 75821-332	UUA-H-HG-C-ATF	21-31"	Task	Glides	Black Self-Skinned Urethane	20" Adjustable	High Profile Aluminum	Chrome	\$918.64
• 75821-334	UUA-H-RC-C-ATF	21-31"	Task	Casters	Black Self-Skinned Urethane	20" Adjustable	High Profile Aluminum	Chrome	\$918.64
• 75821-336	UUA-L-HG-C	16-21"	Task	Glides	Black Self-Skinned Urethane	NA	High Profile Aluminum	Chrome	\$735.90
• 75821-338	UUA-L-RC-C	16-21"	Task	Casters	Black Self-Skinned Urethane	NA	High Profile Aluminum	Chrome	\$735.26
• 75821-340	UUA-M-HG-C	18-25"	Task	Glides	Black Self-Skinned Urethane	NA	High Profile Aluminum	Chrome	\$779.18
• 75821-342	UUA-M-HG-C-ATF	18-25"	Task	Glides	Black Self-Skinned Urethane	20" Adjustable	High Profile Aluminum	Chrome	\$893.40
• 75821-344	UUA-M-RC-C-ATF	18-25"	Task	Casters	Black Self-Skinned Urethane	20" Adjustable	High Profile Aluminum	Chrome	\$892.63
• 75821-346	UUS-H-HG-C	24-31"	Task	Glides	Black Self-Skinned Urethane	21" Affixed	Tubular Steel	Chrome	\$840.57
• 75821-348	UUS-H-RC-C	24-31"	Task	Casters	Black Self-Skinned Urethane	21" Affixed	Tubular Steel	Chrome	\$832.56
• 75821-350	UUS-L-HG-C	17-21"	Task	Glides	Black Self-Skinned Urethane	21" Affixed	Tubular Steel	Chrome	\$798.41
• 75821-352	UUS-M-HG-C	20-27"	Task	Glides	Black Self-Skinned Urethane	21" Affixed	Tubular Steel	Chrome	\$821.25
• 75821-354	UUS-M-RC-C	20-27"	Task	Casters	Black Self-Skinned Urethane	21" Affixed	Tubular Steel	Chrome	\$820.54

*Based on model, Quick Ship products ship next business day (NDS) or in 2 business days (2DS) in limited quantities per order (9 units per NDS order; 15 models per 2DS order).

BIOFIT QUICK SHIP SEATING – WHEN YOU NEED CHAIRS FAST*

FIGHTING INFECTION
THE BIOFIT ADVANTAGE

ETON SERIES ERGONOMIC CHAIR

Available with non-porous PVC chemical-resistant vinyl upholstery which can be cleaned and disinfected.

- Proven design and performance make these chairs a popular choice for today's healthcare and laboratory environments, including ISO 6 clean rooms.
- Features comfortable 18.5" wide x 17" deep x 3" thick seat with four-way contouring and internal seat board bumperguard and streamlined 14.5" wide x 9.5" high backrest.
- Full ergonomic functionality and seat height ranges to accommodate many tasks.



VWR Catalog #	BioFit Part #	Seat Height	Control Style	Foot Type	Upholstery	Footring	Base Style	Metal Finish	Price
• 76404-344	ETSHRTXXXC540	25-32"	Task	Casters	Black Vinyl	21" Affixed	Tubular Steel	Chrome	\$860.93

ETON SERIES ERGONOMIC CHAIR WITH HEIGHT-ADJUSTABLE ARMS

VWR Catalog #	BioFit Part #	Seat Height	Control Style	Foot Type	Upholstery	Footring	Base Style	Metal Finish	Price
• 76404-338	ETALRTXPXC540	17-22"	Task	Casters	Black Vinyl	NA	High Profile Aluminum	Chrome	\$899.40

ETON SERIES ISO 6 CLEAN ROOM CHAIRS

VWR Catalog #	BioFit Part #	Seat Height	Control Style	Foot Type	Upholstery	Footring	Base Style	Metal Finish	Price
• 76404-346	ETSLRTXX6C540	18-22"	Task	Casters	Black Vinyl	21" Affixed	Tubular Steel	Chrome	\$842.91
• 76404-348	ETSMHTXX6C540	21-28"	Task	Glides	Black Vinyl	21" Affixed	Tubular Steel	Chrome	\$863.35

*Based on model, Quick Ship products ship next business day (NDS) or in 2 business days (2DS) in limited quantities per order (9 units per NDS order; 15 models per 2DS order).

BRIDGEPORT SERIES ERGONOMIC CHAIR

Available with non-porous PVC chemical-resistant vinyl upholstery which can be cleaned and disinfected.

- Designed to enhance user comfort and upper body mobility, this series is a flexible seating solution for a variety of healthcare and laboratory workspaces, including ISO 5 clean rooms.
- Concave 18.5" wide x 16.5" deep x 3" thick seat combines with nimble, lumbar-supporting 14.5" wide x 9.5" high backrest to deliver in-demand performance.
- Full ergonomic functionality and seat height ranges to accommodate many tasks.



VWR Catalog #	BioFit Part #	Seat Height	Control Style	Foot Type	Upholstery	Footring	Base Style	Metal Finish	Price
• 76404-310	BTAHRTFXXXC540	22-32"	Task	Casters	Black Vinyl	20" Adjustable	High Profile Aluminum	Chrome	\$888.18
• 76404-318	BTSMRTXXXB540	21-28"	Task	Casters	Black Vinyl	21" Affixed	Tubular Steel	Black Powder Coat	\$665.47
• 76404-328	BTSLRTXXXB540	18-22"	Task	Casters	Black Vinyl	21" Affixed	Tubular Steel	Black Powder Coat	\$633.19

BRIDGEPORT SERIES ISO 5 CLEAN ROOM CHAIRS

VWR Catalog #	BioFit Part #	Seat Height	Control Style	Foot Type	Upholstery	Footring	Base Style	Metal Finish	Price
• 76404-322	BTAHRTFX5CMBK	22-32"	Task	Casters	Black Vinyl	20" Adjustable	High Profile Aluminum	Chrome	\$1,157.26
• 76404-326	BTSHRTXX5CMBK	25-32"	Task	Casters	Black Vinyl	21" Affixed	Tubular Steel	Chrome	\$976.85

**BIOFIT STOOLS: THE SMART CHOICE FOR HEALTHCARE APPLICATIONS –
NOW WITH REDUCED LEAD TIME****

REXFORD RX SERIES STOOL

Available with non-porous PVC chemical-resistant vinyl upholstery which can be cleaned and disinfected.

- 16" diameter x 3.5 thick seat with high-resilience molded foam cushioning and internal seat board bumperguard; 1:10 bleach/water cleanable.
- Fingertip-controlled pneumatic activator ring enables seat height adjustment from any location around the seat perimeter.
- Easy wipe-down chrome finish on metal parts.



VWR Catalog #	BioFit Part #	Seat Height	Control Style	Foot Type	Upholstery	Footring	Base Style	Metal Finish	Price
76404-352	RXSLRSXXXC540	18-22"	Fingertip Height Adjustment Ring	Casters	Black Vinyl	21" Affixed	Tubular Steel	Chrome	\$608.42

*Based on model, Quick Ship products ship next business day (NDS) or in 2 business days (2DS) in limited quantities per order (9 units per NDS order; 15 models per 2DS order).

**Lead time: Ships within 5 business days of order receipt.



Dedicated to Design.
Committed to Quality.

BioFit Engineered Products
PO Box 109
Waterville, OH 43566



1.800.932.5000 | VWR.COM

Prices and product details are current when published and subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2019 Avantor, Inc. All rights reserved.