



VWR® Undercounter and Freestanding/Mobile Glassware Washers

Stainless Steel Construction

Multiple Cycle Selections
and Microprocessor Control

Low Water Requirement



VWR® Undercounter and Freestanding/Mobile Glassware Washers

- Feature stainless steel construction
- Multiple cycle selections and microprocessor control ensure that each load is processed efficiently and with the minimum energy required to meet cleaning requirements
- Low water requirement and precise detergent and neutralizer dosing minimize chemical costs and help protect the environment
- A delayed-start feature allows the washers to operate automatically during off-peak hours
- Undercounter glassware washers fit in a standard cabinet opening and feature adjustable legs
- Freestanding/mobile glassware washers can be installed without occupying cabinet space and are shipped with casters
- Spindle-ready models use an existing deionized water supply and can accommodate either one or two direct injection spindle racks
- Undercounter unit: 61W x 62D x 87–93H cm (24 x 24½ x 34¼–36½")
- Freestanding/mobile unit: 61W x 62D x 97.5H cm (24 x 24½ x 38¾"), minimum height
- Interior dimensions: 55W x 48D x 64H cm (21½ x 19 x 25¼")
- Electrical: 115V, 60Hz, 15A

Description	Cat. No.	Each
Undercounter Glassware Washers		
Washer	82020-922	3496.47
Washer with DI Rinse	82030-586	5668.07
Washer with DI Rinse, Spindle-Ready	82100-002	5842.61
Washer with DI Rinse and DI Pump	97021-360	6429.44
Washer with DI Rinse and DI Pump, Spindle-Ready	97021-364	7038.78
Stainless Steel Base Cabinet for Undercounter Washers	97021-304	1141.66
Freestanding/Mobile Glassware Washers		
Washer	82020-924	4430.96
Washer with DI Rinse	82030-588	6346.36
Washer with DI Rinse, Spindle-Ready	82100-004	6521.68
Washer with DI Rinse and DI Pump	97021-362	7228.59
Washer with DI Rinse and DI Pump, Spindle-Ready	97021-366	7840.07



VWR Lab Planning

Building a new lab? Renovating an existing one? The VWR Furniture team can support you before, during, and after your installation, by offering:

- Comprehensive lab planning assistance
- Choice offerings of quality lab furniture and components
- Complete project supervision from proposal through final installation

Our team of experienced sales and support staff is dedicated to supplying original, innovative solutions - on time and within budget.



Prices, product appearance and specifications are current at the time of printing, subject to change without notice. Availability for certain products may be limited by federal, state, provincial or local licensing requirements. VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in U.S. dollars unless otherwise noted. Offers valid in USA and Canada, void where prohibited by law or company policy, while supplies last. Visit vwr.com to view our privacy policy and additional disclaimers.

VWR, forms of VWR, and the VWR logo and/or design are either registered trademarks ®, trademarks ™, or service marks ™ of VWR International, LLC in the United States and/or other countries. All other marks referenced herein are registered trademarks, trademarks or service marks of their respective owner(s). For a complete list of trademark owners please visit vwr.com.

©2010 VWR International, LLC. All rights reserved.

0210 Lit. No. POD 92715