

avantor™
delivered by VWR™ | VWR.COM

Millipore
Sigma



EASY GRIP FULL CONTROL

Ergonomic HDPE bottles
for acids, bases and solvents

MilliporeSigma is the U.S.
and Canada Life Science
business of Merck KGaA,
Darmstadt, Germany.

Supelco®
Analytical Products

TAKE FULL CONTROL – WITH ONE HAND

New 1 l HDPE bottle and popular
2.5 l HDPE bottle for acids, bases and solvents

Molded for hands. Made for handiness.


Safe and simple handling of hazardous liquids is key to protecting personnel from accidents and health risks. That's why we created the new "easy grip" 1 l HDPE bottle. Just like our popular 2.5 l HDPE bottle, the smaller version has an ergonomically perfected design to ensure a firm grip and full control – no matter how big or small your hands.

Developed with scientists. For scientists.

Today's performance requirements for lab reagents go well beyond product properties. Besides analytical purity, aspects such as handling, safety, ecology and economy play an increasingly crucial role for our customers.

All of these factors are directly influenced by packaging – particularly for acids, bases and solvents. Glass bottles are still the preferred option. As a container for all types of chemicals, glass is a valuable inert material for daily lab use. But there is always the risk of breakage.

At MilliporeSigma, we have been developing the most innovative and practical packaging concepts for many years. Our experts have collaborated with scientists like you to come up with the perfect packaging for solvents, acids and bases: our HDPE (high density polyethylene) bottle range. Developed and used exclusively by MilliporeSigma, it incorporates safety, environmental protection, and cost savings along the entire process chain. Meet our new "easy grip" bottle and feel how it makes your lab work safer, quicker and easier every time.



Benefits you can feel

- Ergonomic shape
- Easy pouring
- Safe handling
- Unbreakable material
- Cost-efficient use

NEW
easy grip
design



PACKED WITH ADDED ADVANTAGES

For maximum user safety and product quality, the bottles are made of a specially treated, high quality HDPE that is extremely durable, inert and shock resistant. What's more, the bottles' unique shapes ensure optimal pouring characteristics and high pressure stability. MilliporeSigma's exclusive S40 thread and lid assure an absolutely tight seal. Moreover, it's designed to provide smooth liquid flow to prevent spills. MilliporeSigma HDPE bottles also protect the environment. Besides cardboard boxes, they require no additional protective material (like polystyrene). So there's much less packaging waste compared to glass bottles.




Scan with
"My M Safety"
Available on
Apple App Store &
Google Play



Ergonomic grip
Optimal handling
and pouring

Smart labels
As with all our products,
the new bottles come
with unique, clear and
complete labels which
include a 2D data matrix
barcode for fast and
easy access to digitalized
product data.

Sharp thread lip
Specially formed for
safe, drip-free pouring



1 l

Height	206 mm
Diameter	101 mm
Minimum volume	1190 ml
Weight	min. 66 g



2.5 l

Height	332 mm
Diameter	125 mm
Minimum volume	2690 ml
Weight	min. 145 g



MilliporeSigma's exclusive S40 thread and lid

Our S40 thread ensures a consistently tight seal thanks to a higher contact pressure. To enhance product safety, every cap is provided with a tamper-evident seal. After the cap is first opened, the tamper-proof ring remains on the bottleneck.

Ergonomic integrated handle

Optimal handling and pouring

Low bottle weight

Easy, safe and economical handling and transportation



For certain chemicals, bottles are colored to protect against UV light.

Special base geometry

Ensures high pressure stability and prevents bulging.

NEW SHAPE. NO CHANGE.

NEW
1 l „easy grip“
HDPE bottle



Previous
1 l HDPE bottle

Our new easy grip 1 l HDPE bottle is made of the same material as its predecessor. So the packaging material is still perfectly compatible with its chemical contents. Also, the new bottle has exactly the same height and diameter, so you won't need to adapt your safety cabinets or shelves.

Ordering information

Acids		
Product	1 l*	2.5 l
Acetic acid 96% for analysis EMSURE®	EM1.00062.1011	EM1.00062.2511
Acetic acid (glacial) 100% anhydrous for analysis EMSURE® ACS, ISO, Reag. Ph Eur	EM1.00063.1011	EM1.00063.2511
Acetic acid (glacial) 100% for analysis EMPARTA® ACS		EM1.01830.2500
Hydrochloric acid 25% for analysis EMSURE®	EM1.00316.1011	EM1.00316.2511
Hydrochloric acid 32% for analysis EMSURE®	EM1.00319.1011	EM1.00319.2511
Hydrochloric acid fuming 37% for analysis EMSURE® ACS, ISO, Reag. Ph Eur	EM1.00317.1011	EM1.00317.2011 ⁺
Hydrofluoric acid 38-40% EMPLURA®	EM1.00337.1000	EM1.00337.2500
Hydrofluoric acid 40% for analysis EMSURE® ISO, Reag. Ph Eur	EM1.00338.1000	EM1.00338.2500
Hydrofluoric acid 48% for analysis EMSURE® ACS, ISO, Reag. Ph Eur	EM1.00334.1000	EM1.00334.2500
Hydrogen peroxide 30% (Perhydrol®) for analysis EMSURE® ISO	EM1.07209.1000	EM1.07209.2500
Hydrogen peroxide 30% (Perhydrol®) (stabilized for higher storage temp.) for analysis EMSURE® ISO	EM1.07210.1000	EM1.07210.2500
ortho-Phosphoric acid 85% for analysis EMSURE® ACS, ISO, Reag. Ph Eur	EM1.00573.1000	EM1.00573.2500
Sulfuric acid 25% for analysis EMSURE®	EM1.00716.1000	
Sulfuric acid 62% for analysis for the determination of fat in cheese (d 1.52) EMSURE®	EM4.80531.1000	EM4.80531.2500
Sulfuric acid 95-97% for analysis EMSURE® ISO	EM1.00731.1011	EM1.00731.2511
Sulfuric acid 96% for the determination of viscosity acc. to DIN EN ISO 307	EM1.08131.1000	EM1.08131.2500

Bases		
Product	1 l*	2.5 l
Ammonia solution 25% for analysis EMSURE® ISO, Reag. Ph Eur	EM1.05432.1011	
Potassium hydroxide solution 47% for analysis EMSURE®	EM1.05545.1000	
Sodium hydroxide solution min. 10% (1.11) for analysis EMSURE®	EM1.05588.1000	
Sodium hydroxide solution min. 27% (1.30) for analysis (for the determination of nitrogen) EMSURE®		EM1.05591.2500
Sodium hydroxide solution about 32% (for the determination of nitrogen) for analysis EMSURE®		EM1.05590.2500
Sodium hydroxide solution about 32% EMPLURA®		EM1.05587.2500
Sodium hydroxide solution min. 45% for analysis EMSURE®		EM1.11360.2500
Sodium hydroxide solution 50% for analysis EMSURE®	EM1.58793.1000	

* New bottle availability is dependent on stock. ⁺ Only available in 2 l

Solvents		
Product	1 l*	2.5 l
Acetone for analysis EMSURE® ACS, ISO, Reag. Ph Eur	EM1.00014.1011	EM1.00014.2511
Acetone for analysis EMPARTA® ACS		EM1.07021.2511
Acetone EMPLURA®	EM8.22251.1000	
tert-Amyl alcohol EMPLURA®	EM8.06193.1000	
1-Butanol EMPLURA®		EM8.22262.2500
2-Butanol EMPLURA®		EM8.22263.2500
Cyclohexane for analysis EMSURE® ACS, ISO, Reag. Ph Eur		EM1.09666.2511
Diethanolamine for analysis EMSURE®	EM1.16205.1000	
N,N-Dimethylformamide for analysis EMSURE® ACS, ISO, Reag. Ph Eur	EM1.03053.1011	EM1.03053.2511
N,N-Dimethylformamide for analysis EMPARTA®	EM1.03034.1011	
N,N-Dimethylformamide EMPLURA®	EM8.22275.1000	EM8.22275.2500
Dimethyl sulfoxide for analysis EMSURE® ACS	EM1.02952.1011	EM1.02952.2511
Ethanol absolute for analysis EMSURE® ACS, ISO, Reag. Ph Eur	EM1.00983.1011	EM1.00983.2511
Ethanol absolute for analysis EMPARTA® ACS		EM1.07017.2511
Ethanol absolute EMPLURA®	EM8.18760.1000	EM8.18760.2500
Ethanol denatured with about 1% methyl ethyl ketone for analysis EMSURE®	EM1.00974.1011	EM1.00974.2511
Ethyl acetate for analysis EMSURE® ACS, ISO, Reag. Ph Eur	EM1.09623.1000	EM1.09623.2511
Ethyl acetate EMPLURA®		EM8.22277.2500
Ethylene glycol for analysis EMSURE® Reag. Ph Eur, Reag. USP	EM1.09621.1000	EM1.09621.2500
Ethylene glycol EMPLURA®	EM1.00949.1000	EM1.00949.2500
Formamide for analysis EMSURE®	EM1.09684.1000	EM1.09984.2500
Formamide EMPLURA®	EM1.04008.1000	EM1.04008.2500
n-Heptane for analysis EMSURE® Reag. Ph Eur		EM1.04379.2511
n-Heptane EMPLURA®		EM1.04365.2511
n-Hexane for analysis EMSURE® ACS		EM1.04367.2511
n-Hexane for analysis EMSURE® ACS, Reag. Ph Eur		EM1.04374.2511
n-Hexane for analysis EMPARTA® ACS		EM1.07023.2511
n-Hexane EMPLURA®		EM1.04368.2511
Isoamyl alcohol (mixture of isomers) EMPLURA®	EM8.22255.1000	EM8.22255.2500
Methanol for analysis EMSURE® ACS, ISO, Reag. Ph Eur	EM1.06009.1011	EM1.06009.2511
Methanol for analysis EMPARTA® ACS		EM1.07018.2511
Methanol EMPLURA®	EM8.22283.1000	EM8.22283.2500
1-Methyl-2-pyrrolidone EMPLURA®	EM8.06072.1000	EM8.06072.2500
1,2-Propanediol EMPLURA®	EM8.22324.1000	
2-Propanol for analysis EMSURE® ACS, ISO, Reag. Ph Eur	EM1.09634.1011	EM1.09634.2511
2-Propanol for analysis EMPARTA® ACS		EM1.07022.2511
2-Propanol EMPLURA®	EM8.18766.1000	EM8.18766.2500
Toluene for analysis EMSURE® ACS, ISO, Reag. Ph Eur		EM1.08325.2511
Toluene for analysis EMPARTA® ACS		EM1.07019.2511

*New bottle availability is dependent on stock.

Safety Accessories	
Product	Ord. No.
Bottle opening key S40/S28	EM1.08801.0001
Safety stand for 2.5 l HDPE bottles	EM9.67213.0001

Full range of different pack sizes
and packaging types available
on VWR.com/MilliporeSigma

Supelco®

Analytical Products

MilliporeSigma
400 Summit Drive
Burlington, MA 01803



We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

© 2022 Merck KGaA, Darmstadt, Germany and its affiliates. All Rights Reserved.
MilliporeSigma, the vibrant M, Chemizorb and Supelco are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners.
Detailed information on trademarks is available via publicly accessible resources.

Lit. No. MS_BR8988EN
08/2022

 **avantor™**
delivered by **VWR™**

VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2023 Avantor, Inc. All rights reserved.

0823 Lit. No.130281W