



Item Number A1009
Item Acetic Acid, Glacial, FCC
CAS Number 64-19-7
Molecular Formula C₂H₄O₂
Molecular Weight 60.05

Assay (C₂H₄O₂)..... 99.5-100.5%
Solidification Point..... Min. 15.6°C
MAXIMUM LIMITS
Lead..... 0.5 mg/kg
Nonvolatile Residue..... 0.005%
Readily Oxidizable Substances..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: **East Coast Facility:**
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843