



# Technical Datasheet

## VWR® SMS Lab Coats - Made in USA

### Intended use

Production, Industrial or Research
------------------------------------

### Characteristics

- Multilayer SMS material protection
- Light fluid and particulate barrier
- Excellent breathability
- Multiple Style Options

SEAM:	Surged Seam Construction
STYLE:	Various, listed in table below
THREAD:	Polyester blend
MATERIAL:	White or Blue SMS 3 layer material 50 gsm
PACKAGED:	Single bag, (poly-lined carton)
QUANTITY:	30 per case
THREAD:	Polyester blend

Size	Cat. No.	Style Description	Case
S	76457-154	white, three pocket, snap, std collar, knit cuffs	30
M	76457-156	white, three pocket, snap, std collar, knit cuffs	30
L	76457-158	white, three pocket, snap, std collar, knit cuffs	30
XL	76457-160	white, three pocket, snap, std collar, knit cuffs	30
2XL	76457-004	white, three pocket, snap, std collar, knit cuffs	30
3XL	76457-006	white, three pocket, snap, std collar, knit cuffs	30
S	76457-008	blue, three pocket, snap, std collar, knit cuffs	30
M	76457-010	blue, three pocket, snap, std collar, knit cuffs	30
L	76457-012	blue, three pocket, snap, std collar, knit cuffs	30
XL	76457-014	blue, three pocket, snap, std collar, knit cuffs	30
2XL	76457-016	blue, three pocket, snap, std collar, knit cuffs	30
3XL	76457-018	blue, three pocket, snap, std collar, knit cuffs	30
4XL	76463-370	blue, three pocket, snap, std collar, knit cuffs	30
5XL	76463-372	blue, three pocket, snap, std collar, knit cuffs	30
6XL	76463-374	blue, three pocket, snap, std collar, knit cuffs	30
7XL	76463-376	blue, three pocket, snap, std collar, knit cuffs	30
S	76457-020	blue, no pocket, snap, std collar open cuffs	30
M	76457-022	blue, no pocket, snap, std collar open cuffs	30
L	76457-024	blue, no pocket, snap, std collar open cuffs	30
XL	76457-026	blue, no pocket, snap, std collar open cuffs	30
2XL	76457-028	blue, no pocket, snap, std collar open cuffs	30
3XL	76457-030	blue, no pocket, snap, std collar open cuffs	30
S	76457-162	white, no pocket, snap, std collar open cuffs	30
M	76457-164	white, no pocket, snap, std collar open cuffs	30



## Technical Datasheet

L	76457-166	white, no pocket, snap, std collar open cuffs	30
XL	76457-168	white, no pocket, snap, std collar open cuffs	30
2XL	76457-032	white, no pocket, snap, std collar open cuffs	30
3XL	76457-034	white, no pocket, snap, std collar open cuffs	30
S	76457-036	white, 3 pocket, snap, knit collar & cuffs	30
M	76457-038	white, 3 pocket, snap, knit collar & cuffs	30
L	76457-040	white, 3 pocket, snap, knit collar & cuffs	30
XL	76457-042	white, 3 pocket, snap, knit collar & cuffs	30
2XL	76457-044	white, 3 pocket, snap, knit collar & cuffs	30
3XL	76457-046	white, 3 pocket, snap, knit collar & cuffs	30
S	76457-048	blue, 3 pocket, snap, knit collar & cuffs	30
M	76457-050	blue, 3 pocket, snap, knit collar & cuffs	30
L	76457-052	blue, 3 pocket, snap, knit collar & cuffs	30
XL	76457-054	blue, 3 pocket, snap, knit collar & cuffs	30
2XL	76457-056	blue, 3 pocket, snap, knit collar & cuffs	30
3XL	76457-058	blue, 3 pocket, snap, knit collar & cuffs	30
S	76457-060	blue, no pockets, standard collar, elastic cuffs	30
M	76457-062	blue, no pockets, standard collar, elastic cuffs	30
L	76457-064	blue, no pockets, standard collar, elastic cuffs	30
XL	76457-066	blue, no pockets, standard collar, elastic cuffs	30
2XL	76457-068	blue, no pockets, standard collar, elastic cuffs	30
3XL	76457-070	blue, no pockets, standard collar, elastic cuffs	30
S	76457-072	white, no pockets, standard collar, elastic cuffs	30
M	76457-074	white, no pockets, standard collar, elastic cuffs	30
L	76457-076	white, no pockets, standard collar, elastic cuffs	30
XL	76457-078	white, no pockets, standard collar, elastic cuffs	30
2XL	76457-080	white, no pockets, standard collar, elastic cuffs	30
3XL	76457-082	white, no pockets, standard collar, elastic cuffs	30

**Packaging:** 30 per Bag /CS 30 PCS Total

**Country of Origin:** USA

Prices, product, and/or services details are current when published and subject to change without notice. Certain products or services may be limited by country, federal, state, provincial, or local regulations. VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. Void where prohibited by law or company policy, while supplies last. Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. Visit [vwr.com](http://vwr.com) to view our privacy policy, trademark owners, and additional disclaimers. © 2020 Avantor, Inc. All rights reserved.

Date of issue: January 2021