



**Item Number** S1764  
**Item** Sucralfate, Fine Powder, USP  
**CAS Number** 54182-58-0  
**Molecular Formula**  $C_{12}H_{54}Al_{16}O_{75}S_8$   
**Molecular Weight** 2086.70

Assay..... 30.0-38.0%  
Aluminum Content..... 15.5-18.5%  
Clarity and Color of Solution..... To pass test

MAXIMUM LIMITS

Acid-neutralizing Capacity..... To pass test  
Chloride..... 0.50%  
Arsenic..... 4 ppm  
Heavy Metals..... 0.002%  
Pyridine and 2-Methylpyridine..... 0.05% (ea.)  
Sucrose Heptasulfate..... To pass test  
Residual solvents..... To pass test

***Spectrum Chemicals & Laboratory Products***

---

**Corporate Headquarters:** 14422 S. San Pedro St.  
Gardena, CA 90248  
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843

**East Coast Facility:** 755 Jersey Ave.  
New Brunswick, NJ 08901