



greiner bio-one

Certificate of Compliance

Product description: Shaker Flask, 1000 ml baffled bottom

Greiner part no.: 679514

Expiration date: Jan. 19, 2028

Greiner lot no.:

1	1	1	2	1	9	-	-
---	---	---	---	---	---	---	---

Physical inspection: passed

Visual inspection: passed

Product has been gamma irradiated and dosimetrically released based upon U.S. Association for the Advancement of Medical Instrumentation (AAMI) Recommended practices. Gamma sterilized to SAL 10^{-6}

The material listed has received the following doses within the precision limits of the dosimetry system employed

Minimum Dosage: 28.5 kGy Maximum Dosage: 35.4 kGy

Isotope Utilized: 60 COBALT-60

Dose Requirement: Min: 18kGy; Max: 30kGy

Pyrogens An extract from the lot contained < 0.005 EU/ml per the Limulus Amoebocyte Lysate (LAL) test method in the FDA *Guidelines on Validation of the Limulus Amoebocyte Lysate Test as an end product Endotoxin test for Human and Animal Parenteral Drugs, Biological Products, and Medical Devices*, December 1987.

RNase-free, DNase-free, DNA-free and Non-pyrogenic.

We hereby certify that this shipment is in accordance with the valid specifications.

23 Apr 2024
Date

Linda Hunter
Quality Assurance

Greiner Bio-One NA Inc
4238 Capital Drive
Monroe, NC 28110, USA