



Item Number SO165
Item Sodium Citrate, Dihydrate, Granular, USP
CAS Number 6132-04-3
Molecular Formula $C_6H_5Na_3O_7 \cdot 2H_2O$
Molecular Weight 294.10

Assay (Anhydrous Basis).....99.0-100.5%
Water..... 10.0-13.0%
Alkalinity..... To pass test

MAXIMUM LIMITS

Tartrate.....To pass test
Heavy Metals.....0.001%
Residual solvents..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: **East Coast Facility:**
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843