


Biohazard and Hazard Bags



The Micronova **Biohazard and Hazard Bags** are constructed from FDA compliant, low density polyethylene. The red and yellow biohazard bags are designed for containment and disposal of biomedical waste with the biohazard symbol and warning message displayed on each biohazard bag. A yellow hazardous material bag is designed for segregation and disposal of other hazardous waste, in areas where red bags are designated biohazard.. Biohazard and Hazard Bags can be used with the BLSSH series of Micronova BinLiner™ stands. Liners gamma irradiated for aseptic areas.



Roll of red biohazard bags with plastic core



- Red or Yellow FDA compliant, low density polyethylene Biohazard and Hazard Bags
- Biohazard symbol in black ink
- 2 mil and 4 mil thicknesses
- Hazardous material bags are designed and labeled for waste segregation that does not go into a biohazard bag
- Smaller bags with a star bottom, ideal for round containers
- Available flat packed or on a plastic core roll
- Good compatibility with common cleaners
- Gamma irradiated bags sterile for aseptic areas

Manufactured for use in Controlled Environments

- ISO 5 thru 9 Cleanrooms
- ISO 9001 Certified Manufacturing Facility
- Good chemical resistance to disinfectants
- Minimum Irradiation Dose: 25 kGy
10⁻⁶ Sterility Assurance Level
- In compliance with AAMI Guidelines for Sterilization (if applicable)

SPECIFICATIONS

Materials of Construction:

Meets FDA Regulation 21 CFR, 177.1520, 3.1
Red or Yellow low density polyethylene

Dimensions & Capacity:

Thickness (mil); Size (in, cm); Capacity (gal, liters)

Biohazard Bags

Red: 2 mil; 25 in x 36 in (63.5 cm x 91.4 cm); 16 gal, 60 liters
2 mil; 34 in x 40 in (86.3 cm x 101.6 cm); 33 gal, 125 liters
4 mil; 40 in x 60 in (101 cm x 152.4 cm); 55 gal, 200 liters

Yellow: 2 mil; 32 in x 36 in (81.3 cm x 91.4 cm); 30 gal, 114 liters
4 mil; 32 in x 40 in (81.3 cm x 101.6 cm); 32 gal, 120 liters

Hazard Bag

Yellow: 2 mil; 24 in x 36 in (61 cm x 91.4 cm); 16 gal, 60 liters



Micronova Manufacturing, Inc.
3431 W. Lomita Boulevard
Torrance, CA 90505-5010 USA

Tel: 310/784-6990
Fax: 310/784-6980

SalesInfo@micronova-mfg.com
www.micronova-mfg.com