

### KELKLAVE Autoclave

#### Autoclave glove for lab applications

The Kelklave Autoclave glove provides heat resistance up to 250°F/121°C.

- Protects up to 250°F/121°C
- White loop-out terry cloth
- Clute cut
- Gauntlett cuff



5" Cuff



11" Cuff

VWR Part #	WL Part #	Cuff
80091-460	422-5	5"
80091-458	422-11	11"

### THERMAL HI-VIS SYNTHETIC SHELL WITH LATEX PALM, A3

#### Freezer glove that protects down to -70°F/-57°C with great dexterity

- Hi-Vis green lightweight synthetic knit shell
- Latex palm provides outstanding grip
- Additional brushed fleece lining material provides heightened thermal properties
- Knit wrist construction prevents dirt and debris from directly contacting the hand
- 13 Gauge
- ANSI Cut Level A3
- ANSI Puncture Level 4
- Available sizes: M-XXL



VWR Part #	WL Part #
10066-726	Y9239TL

### KNIT SHELL WITH FULL NBR COATING AND NITRILE PALM, A4

#### Glass washing glove in lab applications that offers cut and puncture resistance

Coated work gloves provide additional protection in high wear high risk areas of the hand. Y9290 features cut-resistant for additional protection. This is fully dipped with a NBR coating and additionally has a sandy nitrile palm for grip. This double coated glove is an excellent source for gripping and provides extra layers of protection against liquids.

- Stainless steel and polyester construction provides cut and abrasion resistance
- Moderate tactile sensitivity
- Latex free
- ANSI Cut Level A4
- ANSI Puncture Level 4
- Available sizes: S-XXL



VWR Part #	WL Part #
10859-876	Y9290L

### CUT RESISTANT LINER, A4

#### Clean-room cut resistant glove

The cut resistant liner glove is a great option for cut resistance and concealing dirt.

The glove can be used as a liner under another glove for an added protection or worn by itself. An ANSI cut level A4 provides excellent cut resistance, as well as being ambidextrous and latex free. M281 is an ideal glove in industries such as critical environment, clean room, and electronics.



- Stainless steel and polyester construction provides cut and abrasion resistance
- Category III for particulates
- Moderate tactile sensitivity
- Latex free
- Ambidextrous
- ANSI Cut Level A4
- Available sizes: S-XL



VWR Part #	WL Part #
32935-396	M281L

## REUSABLE KNITTED NYLON GLOVE LINER

Wells Lamont Industrial knit liners are made of high-quality, 100% continuous filament. The lightweight glove can be used alone or under a nitrile glove to keep the hand dry from perspiration.

- Stretch nylon construction allows increased flexibility and comfort
- Form fitting nature is perfect for high dexterity applications
- Ambidextrous
- Universal sizing

The **HALF-FINGER** design allows the user tactile in the fingers for movement.

**VWR Part #**      **WL Part #**  
32935-480      M088



The **FULL-FINGER** design provides full protection for those most sensitive glove users.

**VWR Part #**      **WL Part #**  
32935-536      M089



## HIGHLY REUSABLE FULL-FINGER NYLON LINER

### Clean-room glove liner

This highly reusable nylon glove is a great option if you are looking for a glove to withstand up to 30 washings without shrinkage. Made with 100% nylon, this liner offers better mobility and comfort. The M555 is packaged in poly bags and is ideal for critical environment, clean room, electronics, and laboratory areas.

- High performance fiber blended with Lycra® for stretch and reusability
- Category II for particulates
- Good tactile sensitivity
- Ambidextrous
- Available sizes: S-XL

**VWR Part #**      **WL Part #**  
32935-585      M555L



## SPEC-TEC™ STRETCH CUT RESISTANT GLOVE, A2

### Cut resistant glove liner

The stretch cut resistant glove is a great option in providing an added layer of protection performing the job. An ANSI cut level A2 provides good cut resistance, along with being ambidextrous. The glove is launderable up to 12 times without shrinking. M114 is an ideal medical glove for the organ and tissue, pharm and hospital markets.

- High performance fiber blended with Lycra® for stretch and reusability
- Up to 12 laundrings without shrinking
- Good tactile sensitivity
- Ambidextrous
- ANSI Cut Level A2
- Available sizes: S-XL

**VWR Part #**      **WL Part #**  
80091-142      M114L



## SPEC-TEC™ STERILE CUT RESISTANT GLOVE, A1

### Sterile glove liner

This unique thin sterile cut resistant glove is constructed of 100% HPPE fiber that provides comfort underneath a disposable glove. An ANSI cut level 1 provides some cut resistance and is individually wrapped in a sterile package. The M102 is an ideal glove in industries such as medical, organ tissue processing, and pharmaceutical.

- 100% HPPE fiber construction
- Excellent tactile sensitivity
- Sterile
- Latex free
- Full-finger
- Ambidextrous
- ANSI Cut Level A1
- Available sizes: S-XL

**VWR Part #**      **WL Part #**  
80091-122      M102L

