


PDS No. 6501xx	<b>PRODUCT DATA SHEET</b>				Page 1 of 1
Revision 06	96 Well Microplate, PS, U-Bottom				 greiner bio-one
	Greiner Item-No. 6501xx				
Valid for Item-No.:	<b>650101</b>	<b>650160</b> (sterile)	<b>650161</b> (sterile)	<b>650180</b> (sterile)	<b>650185</b> (sterile)

1.	Description / Specification	
1.1	Description	PS Plate, 96 well, clear, solid U-bottom, alphanumeric well coding 650101: standard 650160: physical surface treatment, sterile 650161: standard, sterile 650180: physical surface treatment, with standard lid, sterile 650185: hydrophobic surface (without treatment), with standard lid, sterile
1.2	Dimensions	Plate: length: 127,76 mm (+/- 0,2 mm) width: 85,48 mm (+/- 0,2 mm) curvature: ≤ 200 µm Lid: length: 127,5 mm (+/- 0,2 mm) width: 85,0 mm (+/- 0,2 mm) curvature: +/- 300 µm
1.3	Volume per well	Total volume: 323 µl (mathematically calculated) Working volume: 40 - 280 µl 650160, -180: growth area: 0,35 cm <sup>2</sup>
1.4	Material / Resin	PS (Polystyrene), free of heavy metal
1.5	Colour	Clear
1.6	Sterilization	650101: No 650160, -161, -180, -185: SAL 10 <sup>-3</sup>
1.7	Quality Control	- Raw Material-Control: physical testing - Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	All Item-No.: Free of detectable DNase/RNase, human DNA and pyrogens. 650160, -180, -185: contents non-cytotoxic
2.2	Temperature range	-20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	1000 x g: swinging-bucket rotor
2.5	Chemical Resistance	See homepage: <a href="http://www.gbo.com/bioscience">www.gbo.com/bioscience</a> →Products →Literature →Technical Information→Chemical Resistance of Resins
2.6	Shelf life	650101: N/A 650161, -185: 5 years after month of production 650160, -180: 4 years after month of production
2.7	Other Information	-

3.	Packaging	650101	650160	650161	650180, -185
3.1	Pieces / Bag	10	1	2	1
3.2	Pieces / Box	100	75	100	60
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)			
3.4	Other Information	Certificate of Quality			

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 05	Date 14 July 2014	Date 15 July 2014	Date 15 July 2014	
Date 12.09.11	Name S. Kaelberer	Name Dr. R. Heller	Name A. Schulz	