

## Product Specification

Test	Specification
Assay (NaOH)	49.5 – 52.0 %
Identification	Passes Test
Ammonium Hydroxide Precipitate	≤ 0.010 %
Chloride (Cl)	≤ 0.0025 %
Heavy Metals (as Ag)	≤ 0.0015 %
Iron (Fe)	≤ 0.0005 %
Nitrogen Compounds (as N)	≤ 0.0005 %
Potassium (K)	≤ 0.025 %
Phosphate (PO <sub>4</sub> )	≤ 0.0005 %
Sulfate (SO <sub>4</sub> )	≤ 0.0025 %
Sodium Carbonate (Na <sub>2</sub> CO <sub>3</sub> )	≤ 0.1 %

For Laboratory, Research, or Manufacturing Use  
Appearance (clear, colorless liquid)  
Storage Condition: Protect from air to avoid absorption of carbon dioxide.

Packaging Site: Paris Mfg Ctr & DC