



# Sodium Hydroxide, 10N

## Volumetric Solution

### BAKER ANALYZED Reagent

Product No. 5674  
Specifications current as of:  
Jun 8 2011

TEST	SPECIFICATION
Standardization at 25°C traceable to NIST Standard Reference Material.	
SRM No...	
Appearance (Clear, colorless liquid)	Passes Test
Normality	9.95 - 10.05
Sodium Carbonate ( $\text{Na}_2\text{CO}_3$ )	max. 2.0 %
Trace Impurities (in ppm):	
Chloride (Cl)	max. 5
Heavy Metals (as Pb)	max. 1
Iron (Fe)	max. 0.5
For Laboratory, Research or Manufacturing Use	
IMPORTANT: Material will freeze if stored below 41°F (5°C).	
Protect from air to avoid absorption of carbon dioxide.	