

Calcium Chloride, Pellets
PURIFIED
For Desiccant Use.



Material No.: 1313-01
Revision No.: 0

Product Specification

Test	Specification
Assay (CaCl ₂)	≥ 89 %
Heavy Metals (as Pb)	≤ 0.005 %
Iron (Fe)	≤ 0.02 %
Trace Impurities - Arsenic (As)	≤ 1 ppm

For Laboratory, Research, or Manufacturing Use

Packaging Site: Paris Mfg Ctr & DC


Jamie Ethier
Vice President Global Quality

For questions on this Product Specification please contact Technical Services at 855.282.6867 or +1.610.386.1700
Avantor Performance Materials, LLC
100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone 610.386.1700