



Item Number S1523
Item Sorbitol, Powder, FCC
CAS Number 50-70-4
Molecular Formula C₆H₁₄O₆
Molecular Weight 182.17

Assay (Anhydrous Basis).....91.0-100.5%
pH of a 10% (w/w) solution..... 3.5-7.0

MAXIMUM LIMITS

Lead..... 1 mg/kg
Nickel..... 1 mg/kg
Reducing Sugars.....0.3%
Residue on Ignition..... 0.1%
Water..... 1.5%

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: **East Coast Facility:**
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843