



# MATERIAL SAFETY DATA SHEET

according to the Global Harmonized System

Date of issue: 02/04/2013

Version 1.0

## SECTION 1. Identification

### Product identifier

Product number 820536  
Product name N,N'-Diphenyl-1,4-phenylenediamine for synthesis

### Relevant identified uses of the substance or mixture and uses advised against

Identified uses Chemical for synthesis

### Details of the supplier of the safety data sheet

Company EMD Millipore Corporation | 290 Concord Road, Billerica, MA 01821,  
United States of America | SDS Phone Support: +1-978-715-1335 |  
General Inquiries: +1-978-751-4321 | Monday to Friday, 9:00 AM to  
4:00 PM Eastern Time (GMT-5)  
  
e-mail: mm\_sds@merckgroup.com

Emergency telephone 613-996-6666 CANUTEC (Canada)  
+1-703-527-3887 CHEMTREC (International)  
24 Hours/day; 7 Days/week

## SECTION 2. Hazards identification

### GHS Classification

Skin sensitization, Category 1, H317  
Chronic aquatic toxicity, Category 3, H412  
For the full text of the H-Statements mentioned in this Section, see Section 16.

### GHS-Labeling

*Hazard pictograms*



*Signal Word*  
Warning

*Hazard Statements*

H317 May cause an allergic skin reaction.  
H412 Harmful to aquatic life with long lasting effects.

MATERIAL SAFETY DATA SHEET  
according to the Global Harmonized System

Product number 820536  
Product name N,N'-Diphenyl-1,4-phenylenediamine for synthesis

---

Version 1.0

*Precautionary Statements*

P262 Do not get in eyes, on skin, or on clothing.  
P273 Avoid release to the environment.  
P280 Wear protective gloves.  
P302 + P352 IF ON SKIN: Wash with plenty of soap and water.

**Other hazards**

None known.

---

**SECTION 3. Composition/information on ingredients**

Formula	1,4-(C <sub>6</sub> H <sub>5</sub> NH) <sub>2</sub> C <sub>6</sub> H <sub>4</sub>	C <sub>18</sub> H <sub>16</sub> N <sub>2</sub> (Hill)
CAS-No.	74-31-7	
Molar mass	260.33 g/mol	

**Hazardous ingredients**

*Chemical Name (Concentration)*

CAS-No.

*N,N'-Diphenyl-p-phenylenediamine (>= 90 % - <= 100 % )*  
74-31-7

---

**SECTION 4. First aid measures**

**Description of first-aid measures**

*Inhalation*

After inhalation: fresh air.

*Skin contact*

After skin contact: wash off with plenty of water. Remove contaminated clothing. Consult a physician.

*Eye contact*

After eye contact: rinse out with plenty of water.

*Ingestion*

After swallowing: immediately make victim drink water (two glasses at most). Consult a physician.

Never give anything by mouth to an unconscious person.

**Most important symptoms and effects, both acute and delayed**

Allergic reactions

The following applies to aromatic amines in general: systemic effect: methemoglobinemia with headache, cardiac dysrhythmia, drop in blood pressure, dyspnoea, and spasms, principal symptom: cyanosis (blue discoloration of the blood).

**Indication of any immediate medical attention and special treatment needed**

No information available.

---

**SECTION 5. Fire-fighting measures**

**Extinguishing media**

---

MATERIAL SAFETY DATA SHEET  
according to the Global Harmonized System

Product number 820536  
Product name N,N'-Diphenyl-1,4-phenylenediamine for synthesis

---

Version 1.0

*Suitable extinguishing media*

Water, Carbon dioxide (CO<sub>2</sub>), Foam, Dry powder

*Unsuitable extinguishing media*

For this substance/mixture no limitations of extinguishing agents are given.

**Special hazards arising from the substance or mixture**

Combustible material

Development of hazardous combustion gases or vapors possible in the event of fire.

Fire may cause evolution of:

nitrogen oxides

**Advice for firefighters**

*Special protective equipment for fire-fighters*

Stay in danger area only with self-contained breathing apparatus. Prevent skin contact by keeping a safe distance or by wearing suitable protective clothing.

*Further information*

Suppress (knock down) gases/vapors/mists with a water spray jet. Prevent fire extinguishing water from contaminating surface water or the ground water system.

---

**SECTION 6. Accidental release measures**

**Personal precautions, protective equipment and emergency procedures**

Advice for non-emergency personnel: Avoid substance contact. Avoid inhalation of dusts. Ensure adequate ventilation. Evacuate the danger area, observe emergency procedures, consult an expert.

Advice for emergency responders: Protective equipment see section 8.

**Environmental precautions**

Do not empty into drains.

**Methods and materials for containment and cleaning up**

Cover drains. Collect, bind, and pump off spills.

Observe possible material restrictions (see sections 7 and 10).

Take up dry. Dispose of properly. Clean up affected area. Avoid generation of dusts.

---

**SECTION 7. Handling and storage**

**Precautions for safe handling**

Observe label precautions.

**Conditions for safe storage, including any incompatibilities**

Tightly closed. Dry.

Store at +15°C to +25°C (+59°F to +77°F).

---

**SECTION 8. Exposure controls/personal protection**

**Exposure limit(s)**

Contains no substances with occupational exposure limit values.

**Engineering measures**

MATERIAL SAFETY DATA SHEET  
according to the Global Harmonized System

Product number 820536  
Product name N,N'-Diphenyl-1,4-phenylenediamine for synthesis

---

Version 1.0

Technical measures and appropriate working operations should be given priority over the use of personal protective equipment.

**Individual protection measures**

Protective clothing should be selected specifically for the workplace, depending on concentration and quantity of the hazardous substances handled. The chemical resistance of the protective equipment should be inquired at the respective supplier.

*Hygiene measures*

Immediately change contaminated clothing. Apply skin- protective barrier cream. Wash hands and face after working with substance.

*Eye/face protection*

Safety glasses

*Hand protection*

Chemical-resistant, impervious gloves complying with an approved standard should be worn at all times when handling chemical products if a risk assessment indicates this is necessary.

*Other protective equipment:*

protective clothing

*Respiratory protection*

required when dusts are generated.

Use a properly fitted, air-purifying or air-fed respirator complying with an approved standard if a risk assessment indicates this is necessary. Respirator selection must be based on known or anticipated exposure levels, the hazards of the product and the safe working limits of the selected respirator.

---

**SECTION 9. Physical and chemical properties**

Physical state	crystals
Color	dark gray
Odor	characteristic
Odor Threshold	No information available.
pH	No information available.
Melting point	ca. 135 °C
Boiling point/boiling range	428 - 437 °F (220 - 225 °C) at 1.33 hPa
Flash point	No information available.
Evaporation rate	No information available.
Flammability (solid, gas)	No information available.
Lower explosion limit	No information available.
Upper explosion limit	No information available.

---

MATERIAL SAFETY DATA SHEET  
according to the Global Harmonized System

Product number	820536	Version 1.0
Product name	N,N'-Diphenyl-1,4-phenylenediamine for synthesis	

---

Vapor pressure	No information available.
Relative vapor density	No information available.
Relative density	1.3 g/cm <sup>3</sup> at 68 °F (20 °C)
Water solubility	at 68 °F (20 °C) practically insoluble
Partition coefficient: n-octanol/water	log Pow: 4.04 (calculated) (Lit.) Potential bioaccumulation
Autoignition temperature	No information available.
Decomposition temperature	No information available.
Viscosity, dynamic	No information available.
Explosive properties	No information available.
Bulk density	560 kg/m <sup>3</sup>

---

## SECTION 10. Stability and reactivity

### Reactivity

The following applies in general to flammable organic substances and mixtures: in correspondingly fine distribution, when whirled up a dust explosion potential may generally be assumed.

### Chemical stability

The product is chemically stable under standard ambient conditions (room temperature) .

### Possibility of hazardous reactions

Violent reactions possible with:

Strong oxidizing agents, Strong acids

Caution! In contact with nitrites, nitrates, nitrous acid possible liberation of nitosamines!

### Conditions to avoid

no information available

### Incompatible materials

no information available

### Hazardous decomposition products

in the event of fire: See section 5.

MATERIAL SAFETY DATA SHEET  
according to the Global Harmonized System

Product number 820536  
Product name N,N'-Diphenyl-1,4-phenylenediamine for synthesis

Version 1.0

---

**SECTION 11. Toxicological information**

**Information on toxicological effects**

*Likely route of exposure*

Eye contact, Skin contact, Ingestion

*Acute oral toxicity*

LD50 rat: 2,370 mg/kg (RTECS)

*Acute dermal toxicity*

LD50 rat: 4,500 mg/kg  
(External MSDS)

*Sensitization*

Sensitization test:

Result: positive

(External MSDS)

May cause an allergic skin reaction.

*Genotoxicity in vitro*

Ames test

Salmonella typhimurium

Result: positive

(National Toxicology Program)

Mutagenicity (mammal cell test):

Result: positive

(National Toxicology Program)

Mutagenicity (mammal cell test): chromosome aberration.

Result: negative

(National Toxicology Program)

*Specific target organ systemic toxicity - single exposure*

The substance or mixture is not classified as specific target organ toxicant, single exposure.

*Specific target organ systemic toxicity - repeated exposure*

The substance or mixture is not classified as specific target organ toxicant, repeated exposure.

*Aspiration hazard*

Regarding the available data the classification criteria are not fulfilled.

**Carcinogenicity**

IARC	No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.
OSHA	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.
NTP	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.

MATERIAL SAFETY DATA SHEET  
according to the Global Harmonized System

Product number 820536  
Product name N,N'-Diphenyl-1,4-phenylenediamine for synthesis

---

Version 1.0

ACGIH No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

**Further information**

Under given conditions, contact with nitrites or nitric acid can lead to the formation of nitrosamines, which have shown themselves to be carcinogenic in animal experiments. The following applies to aromatic amines in general: systemic effect: methemoglobinemia with headache, cardiac dysrhythmia, drop in blood pressure, dyspnoea, and spasms, principal symptom: cyanosis (blue discoloration of the blood).

Further data:

Handle in accordance with good industrial hygiene and safety practice.

---

**SECTION 12. Ecological information**

**Ecotoxicity**

*Toxicity to fish*

LC50 *Oryzias latipes* (Orange-red killifish): 36 mg/l; 96 h (External MSDS)

**Persistence and degradability**

*Biodegradability*

0.2 %; 14 d

OECD Test Guideline 301C

Not readily biodegradable.

**Bioaccumulative potential**

*Partition coefficient: n-octanol/water*

log Pow: 4.04

(calculated)

(Lit.) Potential bioaccumulation

**Mobility in soil**

No information available.

**Other adverse effects**

*Additional ecological information*

Discharge into the environment must be avoided.

---

**SECTION 13. Disposal considerations**

The information presented only applies to the material as supplied. The identification based on characteristic(s) or listing may not apply if the material has been used or otherwise contaminated. It is the responsibility of the waste generator to determine the toxicity and physical properties of the material generated to determine the proper waste identification and disposal methods in compliance with applicable regulations. Disposal should be in accordance with applicable regional, national and local laws and regulations.

---

**SECTION 14. Transport information**

**Land transport (DOT)**

Not classified as dangerous in the meaning of transport regulations.

**Air transport (IATA)**

MATERIAL SAFETY DATA SHEET  
according to the Global Harmonized System

Product number 820536  
Product name N,N'-Diphenyl-1,4-phenylenediamine for synthesis

---

Version 1.0

Not classified as dangerous in the meaning of transport regulations.

**Sea transport (IMDG)**

Not classified as dangerous in the meaning of transport regulations.

---

**SECTION 15. Regulatory information**

**United States of America**

**Canada**

**WHMIS Classification**

D2B Toxic Material Causing Other Toxic Effects

Skin sensitizer

This product has been classified according to the hazard criteria of the CPR and the MSDS contains all of the information required by the CPR.

**Notification status**

TSCA: On TSCA Inventory

DSL: All components of this product are on the Canadian DSL.

---

**SECTION 16. Other information**

**Training advice**

Provide adequate information, instruction and training for operators.

**Full text of H-Statements referred to under sections 2 and 3.**

H317 May cause an allergic skin reaction.

H412 Harmful to aquatic life with long lasting effects.

**Key or legend to abbreviations and acronyms used in the safety data sheet**

Used abbreviations and acronyms can be looked up at [www.wikipedia.org](http://www.wikipedia.org).

Date of issue: 02/04/2013

---

The information contained herein is based on the present state of our knowledge. It characterizes the product with regard to appropriate safety precautions. It does not represent a warranty of any product properties and we assume no liability for any loss or injury which may result from the use of this information. Users should conduct their own investigations to determine the suitability of the information.

*All rights reserved. Millipore and the "M" Mark are registered trademarks of Merck KGaA, Darmstadt, Germany.*