

The following list contains the Material Safety Data Sheets you requested. Please scroll down to view the requested MSDS(s).

<u>Product</u>	<u>MSDS</u>	<u>Distributor</u>	<u>Format</u>	<u>Language</u>	<u>Quantity</u>
145201	42432	Hach Company	ROWGHS	English	1
145201	42532	Hach Company	ROWGHS	English	1
145201	42632	Hach Company	ROWGHS	English	1

---

Total Enclosures: 3

World Headquarters  
Hach Company  
P.O.Box 389  
Loveland, CO USA 80539  
(970) 669-3050

MSDS No: M00305

# MATERIAL SAFETY DATA SHEET

---

## 1. CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

**Product Name:** Buffer Solution Hardness 1 pH 10.1 ± 0.1

**Catalog Number:** 42432

Hach Company  
P.O.Box 389  
Loveland, CO USA 80539  
(970) 669-3050

Emergency Telephone Numbers:  
(Medical and Transportation)  
(303) 623-5716 24 Hour Service  
(515)232-2533 8am - 4pm CST

**MSDS Number:** M00305

**Chemical Name:** Not applicable

**CAS Number:** Not applicable

**Additional CAS No. (for hydrated forms):** Not applicable

**Chemical Formula:** Not applicable

**Chemical Family:** Mixture

**Intended Use:** Laboratory Reagent Hardness determination

---

## 2. HAZARDS IDENTIFICATION

**GHS Classification:**

**Hazard categories:** Serious Eye Damage/Eye Irritation: Eye Irrit. 2 Skin Corrosion/Irritation: Skin Irrit. 2 Hazardous to the Aquatic Environment: Aquatic Chronic 3

**GHS Label Elements:**

WARNING



**Hazard statements:** . Causes skin irritation. Causes serious eye irritation. Harmful to aquatic life with long lasting effects.

**Precautionary statements:** Wear protective gloves / protective clothing / eye protection / face protection. If eye irritation persists: Get medical advice/attention. If skin irritation occurs: Get medical advice/attention. Dispose of contents/container according to state, local, federal or national regulations. Handle environmental release according to local, state, federal, provincial requirements. Wear eye protection. IF ON SKIN: Wash with plenty of soap and water. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Call a POISON CENTER or doctor/physician if you feel unwell.

**HMIS:**

**Health:** 1

**Flammability:** 1

**Reactivity:** 0

**Protective Equipment:** X - See protective equipment, Section 8.

**NFPA:**

**Health:** 2

**Flammability:** 1

**Reactivity:** 0

**Symbol:** Not applicable

**WHMIS Hazard Classification:** Class D, Division 2, Subdivision B - Toxic material (other toxic effects)

**WHMIS Symbols:** Other Toxic Effects

---

### 3. COMPOSITION / INFORMATION ON INGREDIENTS

*Hazardous Components according to GHS:*

**Aminomethylpropanol**

*CAS Number:* 124-68-5

*Chemical Formula:* C<sub>4</sub>H<sub>11</sub>NO

*GHS Classification:* Eye Irrit. 2, H319; Skin Irrit. 2, H315; Aquatic Chronic 3, H412

*Percent Range:* < 50

*Percent Range Units:* weight / weight

*PEL:* Not established

*TLV:* Not established

*WHMIS Symbols:* Corrosive/Flammable / Combustible

**Acetic Acid**

*CAS Number:* 64-19-7

*Chemical Formula:* C<sub>2</sub>H<sub>4</sub>O<sub>2</sub>

*GHS Classification:* Flam. Liq 3, H226; Skin Corr. 1A, H314

*Percent Range:* < 10.0

*Percent Range Units:* weight / weight

*PEL:* 10 ppm (25 mg/m<sup>3</sup>)

*TLV:* 10 ppm (25 mg/m<sup>3</sup>)

*WHMIS Symbols:* Corrosive/Flammable / Combustible

*Hazardous Components according to GHS: No*

**Demineralized Water**

*CAS Number:* 7732-18-5

*Chemical Formula:* H<sub>2</sub>O

*GHS Classification:* Non hazardous

*Percent Range:* 30.0 - 40.0

*Percent Range Units:* weight / weight

*PEL:* Not established

*TLV:* Not established

*WHMIS Symbols:* Not applicable

**Magnesium Disodium EDTA Solution**

*CAS Number:*

*Chemical Formula:*

*GHS Classification:* Not applicable

*Percent Range:* < 0.5

*Percent Range Units:* weight / weight

*PEL:* Not established

*TLV:* Not established

*WHMIS Symbols:* Other Toxic Effects

---

### 4. FIRST AID MEASURES

*General Information:* In the event of exposure, show this Material Safety Data Sheet and label (where possible) to a doctor.

*Advice to doctor:* Treat symptomatically.

*Eye Contact:* Immediately flush eyes with water for 15 minutes. Call physician if irritation develops.

*Skin Contact (First Aid):* Wash skin with plenty of water. Call physician if irritation develops.

*Inhalation:* Remove to fresh air. Give artificial respiration if necessary.

*Ingestion (First Aid):* Give large quantities of water. If you feel unwell, contact a physician.

---

### 5. FIRE FIGHTING MEASURES

**Flammable Properties:** Material is not classified as flammable according to GHS criteria. Can burn in fire, releasing toxic vapors.

**Fire Fighting Instruction:** As in any fire, wear self-contained breathing apparatus pressure-demand and full protective gear.

**Extinguishing Media:** Water. Dry chemical. Carbon dioxide. Alcohol foam.

**Extinguishing Media NOT To Be Used:** Not applicable

**Fire / Explosion Hazards:** May react violently with: strong oxidizers

**Hazardous Combustion Products:** Toxic fumes of: nitrogen oxides. carbon monoxide, carbon dioxide.

---

## 6. ACCIDENTAL RELEASE MEASURES

### **Spill Response Notice:**

Only persons properly qualified to respond to an emergency involving hazardous substances may respond to a spill according to federal regulations (OSHA 29 CFR 1910.120(a)(v)) and per your company's emergency response plan and guidelines/procedures. See Section 13, Special Instructions for disposal assistance. Outside of the US, only persons properly qualified according to state or local regulations should respond to a spill involving chemicals.

**Containment Technique:** Absorb spilled liquid with non-reactive sorbent material. Dike large spills to keep spilled material from entering sewage and drainage systems or bodies of water.

**Clean-up Technique:** If permitted by regulation, Cover spilled material with a dry acid, such as citric or boric. Scoop up slurry into a large beaker. Dilute with a large excess of water. Adjust to a pH between 6 and 9 with an acid, such as sulfuric or citric. Flush reacted material to the drain with a large excess of water. Decontaminate the area of the spill with a weak acid solution. Otherwise, Pick up spill for disposal and place in a closed container. Dispose of in accordance with local, state and federal regulations or laws.

**Evacuation Procedure:** Evacuate as needed to perform spill clean-up. If conditions warrant, increase the size of the evacuation.

**DOT Emergency Response Guide Number:** Not applicable

---

## 7. HANDLING AND STORAGE

**Handling:** Avoid contact with eyes skin clothing. Do not breathe mist or vapors. Wash thoroughly after handling. Maintain general industrial hygiene practices when using this product.

**Storage:** Store away from: oxidizers. Protect from: heat

**Flammability Class:** Class IIIB

---

## 8. EXPOSURE CONTROLS / PERSONAL PROTECTION

**Engineering Controls:** Maintain general industrial hygiene practices when using this product.

### **Personal Protective Equipment:**

**Eye Protection:** chemical splash goggles

**Skin Protection:** lab coat nitrile gloves. In the EU, the selected gloves must satisfy the specifications of EU Directive 89/686/EEC and standard EN 374 derived from it.

**Inhalation Protection:** adequate ventilation

**Precautionary Measures:** Avoid contact with: eyes skin clothing. Do not breathe: mist/vapor. Wash thoroughly after handling. Keep away from: oxidizers

**TLV:** Not established

**PEL:** Not established

**For Occupational Exposure Limits (OEL) for ingredients, see section 3 - Composition/Information on Ingredients.:**

---

## 9. PHYSICAL AND CHEMICAL PROPERTIES

**Appearance:** Clear, yellow liquid

**Physical State:** Liquid

**Molecular Weight:** Not applicable

**Odor:** Vinegar

**Odor Threshold:** Not established

**pH:** 10.0

### **Metal Corrosivity:**

**Corrosivity Classification:** Not classified as corrosive to metals according to GHS criteria.

**Steel:** 0.002 in/yr

**Aluminum:** Not determined

**Specific Gravity/ Relative Density (water = 1; air =1):** 1.033

**Viscosity:** Not determined

**Solubility:**

**Water:** Soluble

**Acid:** Soluble

**Other:** Not determined

**Partition Coefficient (n-octanol / water):** Not applicable

**Coefficient of Water / Oil:** Not determined

**Melting Point:** -16 °C (3 °F)

**Decomposition Temperature:** Not determined

**Boiling Point:** 104.5 °C (220 °F)

**Vapor Pressure:** 23 mm Hg @ 25 °C (77 °F)

**Vapor Density (air = 1):** 0.6

**Evaporation Rate (water = 1):** 0.97

**Volatile Organic Compounds Content:** Not determined

**Flammable Properties:** Material is not classified as flammable according to GHS criteria. Can burn in fire, releasing toxic vapors.

**Flash Point:** > 97.2 °C (> 207 °F)

**Method:** Closed cup

**Flammability Limits:**

**Lower Explosion Limits:** Not determined

**Upper Explosion Limits:** Not determined

**Autoignition Temperature:** Not determined

**Explosive Properties:**

Not classified according to GHS criteria.

**Oxidizing Properties:**

Not classified according to GHS criteria.

**Reactivity Properties:**

Not classified as self-reactive, pyrophoric, self-heating or emitting flammable gases in contact with water according to GHS criteria.

**Gas under Pressure:**

Not classified according to GHS criteria.

Not applicable

---

## 10. STABILITY AND REACTIVITY

**Chemical Stability:** Stable when stored under proper conditions.

**Mechanical Impact:** None reported

**Static Discharge:** None reported.

**Reactivity / Incompatibility:** May react violently in contact with: oxidizers

**Hazardous Decomposition:** Toxic fumes of: nitrogen oxides carbon dioxide carbon monoxide

**Conditions to Avoid:** Extreme temperatures

---

## 11. TOXICOLOGICAL INFORMATION

**Toxicokinetics, Metabolism and Distribution:**

Acetic acid exists only in a ionized form, therefore it is not included in determining the hazards of the mixture.

**Toxicologically Synergistic Products:** None reported

**Acute Toxicity:** Acute Toxicity Estimate (ATE) - Calculated from Ingredient Toxicity Data Route Data Given Below

Based on classification principles, the classification criteria are not met.

Oral Rat LD50 = 5566 mg/kg

Dermal Rabbit LD50 = 3839 mg/kg

**Specific Target Organ Toxicity - Single Exposure (STOT-SE):** Based on classification principles, the classification criteria are not met.

**Specific Target Organ Toxicity - Repeat Exposure (STOT-RE):** Based on classification principles, the classification criteria are not met.

**Skin Corrosion/Irritation:** Irritating to skin.

**Eye Damage:** Irritating to eyes.

**Sensitization:** Based on classification principles, the classification criteria are not met.

**CMR Effects/Properties (carcinogenic, mutagenic or toxic to reproduction):** Based on classification principles, the classification criteria are not met.

This product does NOT contain any IARC listed chemicals.  
This product does NOT contain any NTP listed chemicals.  
This product does NOT contain any OSHA listed carcinogens.

**Symptoms/Effects:**

**Ingestion:** May cause: abdominal pain gastrointestinal tract irritation nausea diarrhea vomiting

**Inhalation:** No effects anticipated Large doses may cause: respiratory tract irritation

**Skin Absorption:** May be harmful if absorbed through skin.

**Chronic Effects:** Chronic overexposure may cause eye irritation skin irritation

**Medical Conditions Aggravated:** Pre-existing: Eye conditions Skin conditions Respiratory conditions

---

## 12. ECOLOGICAL INFORMATION

**Product Ecological Information:** --

Do not place in landfill. Recycle appropriately. Do not release into the environment.

Method Used for Estimation of Aquatic Toxicity of Mixture CEPA Categorization: Not Persistent or Bioaccumulative.

Not inherently toxic to aquatic organisms. M-factor (Multiplier) for highly toxic ingredients: 1

**Ingredient Ecological Information:** Aminomethylpropanol: 96 hr Pleuronectes platessa LC50 = 184 mg/L; 48 hr Daphnia magna EC50 = 193 mg/L; 24 hr Daphnia magna EC50 = 65 mg/L; 72 hr Scenedesmus subspicatus EC50 = 520 mg/L

CEPA categorization for ingredients are as follows:

Aminomethylpropanol: Not persistent, bioaccumulative or inherently toxic to aquatic organisms.

---

## 13. DISPOSAL CONSIDERATIONS

**EPA Waste ID Number:** Not applicable

**Special Instructions (Disposal):** Dilute to 3 to 5 times the volume with cold water. Adjust to a pH between 6 and 9 with an acid, such as sulfuric or citric. If permitted by regulation, Open cold water tap completely, slowly pour the reacted material to the drain. Allow cold water to run for 5 minutes to completely flush the system. Otherwise, Check with national, local municipal and state authorities and waste contractors for pertinent local information on the disposal of this article.

**Empty Containers:** Rinse three times with an appropriate solvent. Collect rinsate and dispose of according to local, state or federal regulations. In the US, rinsate from empty containers is classified as hazardous waste and should be disposed of at an E.P. A. approved facility. Rinsate from empty containers may contain sufficient product to require disposal as hazardous waste.

**NOTICE (Disposal):** These disposal guidelines are based on federal regulations and may be superseded by more stringent state or local requirements. Please consult your local environmental regulators for more information. In Europe: Chemical and analysis solutions must be disposed of in compliance with the respective national regulations. Product packaging must be disposed of in compliance with the country-specific regulations or must be passed to a packaging return system.

---

## 14. TRANSPORT INFORMATION

**D.O.T.:**

**D.O.T. Proper Shipping Name:** Not Currently Regulated

--

**Hazard Class:** NA

**Subsidiary Risk:** NA

**ID Number:** NA

**Packing Group:** NA

**T.D.G.:**

**Proper Shipping Name:** Not Currently Regulated

--

**Hazard Class:** NA

**Subsidiary Risk:** NA

**UN Number/PIN:** NA

**Packing Group:** NA

**I.C.A.O.:**

**I.C.A.O. Proper Shipping Name:** Not Currently Regulated

--

**Hazard Class:** NA

**Subsidiary Risk:** NA

**ID Number:** NA

**Packing Group:** NA

**I.M.O.:**

**Proper Shipping Name:** Not Currently Regulated

--

**Hazard Class:** NA

**Subsidiary Risk:** NA

**ID Number:** NA

**Packing Group:** NA

**Additional Information:** There is a possibility that this product could be contained in a reagent set or kit composed of various compatible dangerous goods. If the item is NOT in a set or kit, the classification given above applies. If the item IS part of a set or kit, the classification would change to the following: UN3316 Chemical Kit, Class 9, PG II or III. If the item is not regulated, the Chemical Kit classification does not apply.

---

## 15. REGULATORY INFORMATION

### **U.S. Federal Regulations:**

**O.S.H.A.:** This product meets the criteria for a hazardous substance as defined in the Hazard Communication Standard. (29 CFR 1910.1200)

### **E.P.A.:**

**S.A.R.A. Title III Section 311/312 Categorization (40 CFR 370):** Immediate (Acute) Health Hazard

**S.A.R.A. Title III Section 313 (40 CFR 372):** This product does NOT contain any chemical subject to the reporting requirements of Section 313 of Title III of SARA.

--

**302 (EHS) TPQ (40 CFR 355):** Not applicable

**304 CERCLA RQ (40 CFR 302.4):** Acetic acid 5000 lbs.

**304 EHS RQ (40 CFR 355):** Not applicable

**Clean Water Act (40 CFR 116.4):** Acetic acid - RQ 5000 lbs.

**RCRA:** Contains no RCRA regulated substances.

### **State Regulations:**

**California Prop. 65:** No Prop. 65 listed chemicals are present in this product.

**Identification of Prop. 65 Ingredient(s):** Not applicable

**California Perchlorate Rule CCR Title 22 Chap 33:** Not applicable

**Trade Secret Registry:** Not applicable

### **National Inventories:**

**U.S. Inventory Status:** All ingredients in this product are listed on the TSCA 8(b) Inventory (40 CFR 710).

**CAS Number:** Not applicable

**Canadian Inventory Status:** All ingredients of this product are DSL/NDSL Listed.

**EEC Inventory Status:** All ingredients used to make this product are listed on EINECS / ELINCS.

**Australian Inventory (AICS) Status:** All ingredients are listed.

**New Zealand Inventory (NZIoC) Status:** All components either listed or exempt.

**Korean Inventory (KECI) Status:** All components of this product are either listed, listed as the anhydrous compound or exempt.

**Japan (ENCS) Inventory Status:** All components either listed or exempt.

**China (PRC) Inventory (MEP) Status:** All components either listed or exempt.

---

## 16. OTHER INFORMATION

**References:** 29 CFR 1900 - 1910 (Code of Federal Regulations - Labor). Air Contaminants, Federal Register, Vol. 54, No. 12. Thursday, January 19, 1989. pp. 2332-2983. Fire Protection Guide on Hazardous Materials, 10th Ed. Quincy, MA: National Fire Protection Fire Protection Guide on Hazardous Materials, 10th Ed. Quincy, MA: National Fire Protection Association, 1991. In-house information. List of Dangerous Substances Classified in Annex I of the EEC Directive (67/548) - Classification, Packaging and Labeling of Dangerous Substances, Amended July 1992. Technical Judgment. TLV's Threshold Limit Values and Biological Exposure Indices for 1992-1993. American Conference of Governmental Industrial Hygienists, 1992.

**Complete Text of H phrases referred to in Section 3:** H315 Causes skin irritation. H319 Causes serious eye irritation. H412 Harmful to aquatic life with long lasting effects.

**Revision Summary:** . Substantial revision to comply with EU Reg 1272/2008, Reg 1907/2006 and UN GHS (ST/SG/AC.10/36/Add.3).

### **Date of MSDS Preparation:**

**Day:** 11

**Month:** December

**Year:** 2013

**MSDS Prepared:** MSDS prepared by Product Compliance Department extension 3350

**CCOHS Evaluation Note:** It is offered under exemption from WHMIS labeling as specified in the Controlled Products Regulation (CPR) Section 17. It is offered under the interim policy that was established by Health Canada permitting use of GHS-formatted safety data sheets in Canada prior to revision of CPR to GHS. This product has been classified and labeled in accordance with the requirements of GHS (ST/SG/AC.10/36/Add.3). This SDS has been prepared in accordance with the requirements of GHS (ST/SG/AC.10/36/Add.3).

---

**Legend:**

NA - Not Applicable	w/w - weight/weight
ND - Not Determined	w/v - weight/volume
NV - Not Available	v/v - volume/volume

**USER RESPONSIBILITY:** Each user should read and understand this information and incorporate it in individual site safety programs in accordance with applicable hazard communication standards and regulations.

**THE INFORMATION CONTAINED HEREIN IS BASED ON DATA CONSIDERED TO BE ACCURATE. HOWEVER, NO WARRANTY IS EXPRESSED OR IMPLIED REGARDING THE ACCURACY OF THESE DATA OR THE RESULTS TO BE OBTAINED FROM THE USE THEREOF.**

**HACH COMPANY ©2014**



World Headquarters  
Hach Company  
P.O.Box 389  
Loveland, CO USA 80539  
(970) 669-3050

MSDS No: M00635

# MATERIAL SAFETY DATA SHEET

## 1. CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

**Product Name:** ManVer ® Hardness Indicator  
**Catalog Number:** 42532

Hach Company  
P.O.Box 389  
Loveland, CO USA 80539  
(970) 669-3050

Emergency Telephone Numbers:  
(Medical and Transportation)  
(303) 623-5716 24 Hour Service  
(515)232-2533 8am - 4pm CST

**MSDS Number:** M00635  
**Chemical Name:** Not applicable  
**CAS Number:** Not applicable  
**Additional CAS No. (for hydrated forms):** Not applicable  
**Chemical Formula:** Not applicable  
**Chemical Family:** Not applicable  
**Intended Use:** Laboratory Reagent Indicator for hardness

## 2. HAZARDS IDENTIFICATION

**GHS Classification:**

**Hazard categories:** Flammable Liquids: Flam. Liq. 3 Corrosive to Metals: Met. Corr. 1 Carcinogenicity: Carc. 2  
Respiratory or Skin Sensitization: Skin Sens.1

**GHS Label Elements:**

WARNING



**Hazard statements:** Suspected of causing cancer. Flammable liquid and vapour. May be corrosive to metals. May cause an allergic skin reaction. Suspected of causing cancer by inhalation.

**Precautionary statements:** Keep only in original container. Use personal protective equipment as required. IF ON SKIN: Wash with plenty of soap and water. IF exposed or concerned: Get medical advice/attention. If skin irritation or rash occurs: Get medical advice/attention. Call a POISON CENTER or doctor/physician if you feel unwell. Keep away from heat/sparks/open flames/hot surfaces. - No smoking. Avoid breathing dust/fume/gas/mist/vapours/spray. IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.

**HMIS:**

**Health:** 2\*

**Flammability:** 3

**Reactivity:** 0

**Protective Equipment:** X - See protective equipment, Section 8.

**NFPA:**

**Health:** 1

**Flammability:** 3

**Reactivity:** 0

**Symbol:** Not applicable

**WHMIS Hazard Classification:** Class D, Division 1, Subdivision B - Toxic material (immediate effects) Class B, Division 2 - Flammable liquids

**WHMIS Symbols:** Acute Poison Flammable / Combustible

---

### 3. COMPOSITION / INFORMATION ON INGREDIENTS

*Hazardous Components according to GHS:*

**Isopropanol**

*CAS Number:* 67-63-0

*Chemical Formula:* C<sub>3</sub>H<sub>8</sub>O

*GHS Classification:*

*Percent Range:* < 5.0

*Percent Range Units:* volume / volume

*PEL:* 400 ppm (980 mg/m<sup>3</sup>)

*TLV:* 200 ppm (492 mg/m<sup>3</sup>)

*WHMIS Symbols:* Flammable / Combustible Other Toxic Effects

**Hydroxylamine Hydrochloride**

*CAS Number:* 5470-11-1

*Chemical Formula:* NH<sub>2</sub>OH·HCl

*GHS Classification:* Met. Corr. 1, H290; Carc. 2, H351; Acute Tox. 4, H312; Acute Tox. 3-Orl, H301; STOT RE 2, H373; Eye Irrit. 2, H319; Skin Irrit. 2, H315; Skin Sens. 1, H317; Aquatic Acute 1, H400

*Percent Range:* 1.0 - 10.0

*Percent Range Units:* weight / volume

*PEL:* 15 mg/m<sup>3</sup> as inhalable dust; 5 mg/m<sup>3</sup> as respirable dust

*TLV:* 10 mg/m<sup>3</sup> as inhalable dust; 3 mg/m<sup>3</sup> as respirable dust

*WHMIS Symbols:* Acute Poison

**Calmagite**

*CAS Number:* 3147-14-6

*Chemical Formula:* C<sub>17</sub>H<sub>14</sub>N<sub>2</sub>O<sub>5</sub>S·H<sub>2</sub>O

*GHS Classification:* Skin Irrit 2, H315; Eye Irrit 2A, H319; STOT Single 3, H335

*Percent Range:* < 0.1

*Percent Range Units:* weight / volume

*PEL:* Not established

*TLV:* Not established

*WHMIS Symbols:* Not applicable

*Hazardous Components according to GHS:* No

**Propylene Glycol**

*CAS Number:* 57-55-6

*Chemical Formula:* C<sub>3</sub>H<sub>8</sub>O<sub>2</sub>

*GHS Classification:* Not applicable

*Percent Range:* 90.0 - 100.0

*Percent Range Units:* volume / volume

*PEL:* Not established

*TLV:* Not established

*WHMIS Symbols:* Not applicable

---

### 4. FIRST AID MEASURES

*General Information:* In the event of exposure, show this Material Safety Data Sheet and label (where possible) to a doctor.

*Advice to doctor:* Treat symptomatically.

*Eye Contact:* Immediately flush eyes with water for 15 minutes. Call physician.

*Skin Contact (First Aid):* Wash skin with soap and plenty of water. Call physician if irritation develops.

*Inhalation:* Remove to fresh air.

*Ingestion (First Aid):* Do not induce vomiting. Give 1-2 glasses of water. Call physician immediately. Never give anything by mouth to an unconscious person.

---

## 5. FIRE FIGHTING MEASURES

**Flammable Properties:** Combustion generates toxic fumes.

**Fire Fighting Instruction:** As in any fire, wear self-contained breathing apparatus pressure-demand and full protective gear.

**Extinguishing Media:** Alcohol foam.

**Extinguishing Media NOT To Be Used:** Not applicable

**Fire / Explosion Hazards:** May react violently with: strong oxidizers Do not expose to sparks or other ignition sources.

**Hazardous Combustion Products:** Toxic fumes of: chlorides carbon monoxide, carbon dioxide.

---

## 6. ACCIDENTAL RELEASE MEASURES

**Spill Response Notice:**

Only persons properly qualified to respond to an emergency involving hazardous substances may respond to a spill according to federal regulations (OSHA 29 CFR 1910.120(a)(v)) and per your company's emergency response plan and guidelines/procedures. See Section 13, Special Instructions for disposal assistance. Outside of the US, only persons properly qualified according to state or local regulations should respond to a spill involving chemicals.

**Containment Technique:** Releases of this material may contaminate the environment. Remove all combustible material from spill area. Remove all ignition and spark-creating sources from the spill area. Cover spilled liquid with a commercially available flammable liquid sorbent such as vapor barrier blanket or activated carbon to avoid evolution of fumes. Vapors may travel to a source of ignition and flash back. May be ignited by: heat, sparks, or flames. Dike the material to create a barrier to combustibles.

**Clean-up Technique:** Eliminate all sources of ignition. Do not breathe the fumes. Use only non-sparking tools. Cover spilled material with an alkali, such as soda ash or sodium bicarbonate. If permitted by regulation, Scoop up slurry into a large beaker. Dilute with a large excess of water. Filter to remove solids. Flush the spilled material to the drain with a large excess of water. Decontaminate the area of the spill with a soap solution. Otherwise, Dispose of in accordance with local, state and federal regulations or laws.

**Evacuation Procedure:** Evacuate general area (50 foot radius or as directed by your facility's emergency response plan) when: any quantity is spilled. If conditions warrant, increase the size of the evacuation.

**DOT Emergency Response Guide Number:** 132

---

## 7. HANDLING AND STORAGE

**Handling:** Avoid contact with eyes skin clothing Do not breathe mist or vapors. Wash thoroughly after handling. Maintain general industrial hygiene practices when using this product.

**Storage:** Keep away from: oxidizers Protect from: sparks, flames and other ignition sources

**Flammability Class:** Class IC

---

## 8. EXPOSURE CONTROLS / PERSONAL PROTECTION

**Engineering Controls:** Maintain general industrial hygiene practices when using this product.

**Personal Protective Equipment:**

**Eye Protection:** safety glasses with top and side shields

**Skin Protection:** disposable latex gloves In the EU, the selected gloves must satisfy the specifications of EU Directive 89/686/EEC and standard EN 374 derived from it. lab coat

**Inhalation Protection:** adequate ventilation

**Precautionary Measures:** Avoid contact with: eyes skin clothing Do not breathe: mist/vapor Wash thoroughly after handling. Keep away from: oxidizers Protect from: sparks, flames and other ignition sources

**TLV:** Not established

**PEL:** Not established

**For Occupational Exposure Limits (OEL) for ingredients, see section 3 - Composition/Information on Ingredients.:**

---

## 9. PHYSICAL AND CHEMICAL PROPERTIES

**Appearance:** Dark red liquid

**Physical State:** Liquid

**Molecular Weight:** Not applicable

**Odor:** Fruity

**Odor Threshold:** Not available

**pH:** 1.09

**Metal Corrosivity:**

**Corrosivity Classification:** Classified as corrosive to metals.

**Steel:** 0.288 in/yr

**Aluminum:** 0.001 in/yr

**Specific Gravity/ Relative Density (water = 1; air =1):** 1.01

**Viscosity:** Not determined

**Solubility:**

**Water:** Soluble

**Acid:** Soluble

**Other:** Not determined

**Partition Coefficient (n-octanol / water):** Not applicable

**Coefficient of Water / Oil:** Not applicable

**Melting Point:** Not determined

**Decomposition Temperature:** Not determined

**Boiling Point:** 118°C

**Vapor Pressure:** Not determined

**Vapor Density (air = 1):** Not determined

**Evaporation Rate (water = 1):** 0.05

**Volatile Organic Compounds Content:** Not determined

**Flammable Properties:** Combustion generates toxic fumes.

**Flash Point:** 25.7°C (78.3°F)

**Method:** Closed cup

**Flammability Limits:**

**Lower Explosion Limits:** Not determined

**Upper Explosion Limits:** Not determined

**Autoignition Temperature:** Not determined

**Explosive Properties:**

Not classified according to GHS criteria.

**Oxidizing Properties:**

Not classified according to GHS criteria.

**Reactivity Properties:**

Not classified as self-reactive, pyrophoric, self-heating or emitting flammable gases in contact with water according to GHS criteria.

**Gas under Pressure:**

Not classified according to GHS criteria.

---

## 10. STABILITY AND REACTIVITY

**Chemical Stability:** Stable when stored under proper conditions.

**Mechanical Impact:** None reported

**Static Discharge:** None reported.

**Reactivity / Incompatibility:** Incompatible with: oxidizers

**Hazardous Decomposition:** Toxic fumes of: chlorides carbon monoxide carbon dioxide

**Conditions to Avoid:** Contact with heat, sparks, open flames or other ignition sources. Heating to decomposition.

---

## 11. TOXICOLOGICAL INFORMATION

**Toxicokinetics, Metabolism and Distribution:** No information available for mixture.

**Toxicologically Synergistic Products:** None reported

**Acute Toxicity:** Acute Toxicity Estimate (ATE) - Calculated from Ingredient Toxicity Data

ATE Oral Rat LD50 = 2701 mg/kg

ATE Dermal Rat LD50 = 20921 mg/kg

**Specific Target Organ Toxicity - Single Exposure (STOT-SE):** Based on classification principles, the classification criteria are not met.

**Specific Target Organ Toxicity - Repeat Exposure (STOT-RE):** Based on classification principles, the classification criteria are not met.

**Skin Corrosion/Irritation:** Irritating to skin.

**Eye Damage:** Based on classification principles, the classification criteria are not met.

**Sensitization:** Skin Sensitizer

**CMR Effects/Properties (carcinogenic, mutagenic or toxic to reproduction):** Contains Listed Carcinogen

Contains < 5% Hydroxylamine Hydrochloride  
This product does NOT contain any IARC listed chemicals.  
This product does NOT contain any NTP listed chemicals.  
This product does NOT contain any OSHA listed carcinogens.

**Symptoms/Effects:**

**Ingestion:** Very large doses may cause: central nervous system effects dizziness drowsiness incoordination headache abdominal cramps rapid pulse and respirations convulsions Hydroxylamine Hydrochloride causes a decreased supply of oxygen to the tissues, blue discoloration of the skin, convulsions, drop in blood pressure and coma.

**Inhalation:** May cause: irritation of nose and throat

**Skin Absorption:** Will be absorbed through the skin. Effects similar to those of ingestion

**Chronic Effects:** Chronic overexposure may cause damage to red blood cells

**Medical Conditions Aggravated:** Pre-existing: Eye conditions Skin conditions Respiratory conditions

---

## 12. ECOLOGICAL INFORMATION

**Product Ecological Information:** --

No ecological data available for this product. Mobility in soil: No data available

**Ingredient Ecological Information:** Hydroxylamine Hydrochloride: LC50 Leuciscus idus 48h = 1-10 mg/l;

CEPA Statement: Hydroxylamine Hydrochloride, Isopropanol and Calmagite: Persistent, not bioaccumulative or inherently toxic to aquatic organisms; Propylene Glycol: Not persistent, not bioaccumulative or inherently toxic to aquatic organisms.

---

## 13. DISPOSAL CONSIDERATIONS

**EPA Waste ID Number:** D001 D002

**Special Instructions (Disposal):** Dilute to 3 to 5 times the volume with cold water. Adjust to a pH between 6 and 9 with an alkali, such as soda ash or sodium bicarbonate. Open cold water tap completely, slowly pour the reacted material to the drain.

**Empty Containers:** Rinse three times with an appropriate solvent. Dispose of empty container as normal trash.

**NOTICE (Disposal):** These disposal guidelines are based on federal regulations and may be superseded by more stringent state or local requirements. Please consult your local environmental regulators for more information. In Europe: Chemical and analysis solutions must be disposed of in compliance with the respective national regulations. Product packaging must be disposed of in compliance with the country-specific regulations or must be passed to a packaging return system.

---

## 14. TRANSPORT INFORMATION

**D.O.T.:**

**D.O.T. Proper Shipping Name:** Flammable Liquid, Corrosive, N.O.S.

(<10% Isopropanol/Hydroxylamine Hydrochloride Solution)

**Hazard Class:** 3

**Subsidiary Risk:** 8

**ID Number:** UN2924

**Packing Group:** III

**T.D.G.:**

**Proper Shipping Name:** Flammable Liquid, Corrosive, N.O.S.

(<10% Isopropanol/Hydroxylamine Hydrochloride Solution)

**Hazard Class:** 3.3

**Subsidiary Risk:** 8

**UN Number/PIN:** 2924

**Packing Group:** III

**I.C.A.O.:**

**I.C.A.O. Proper Shipping Name:** FLAMMABLE LIQUID, CORROSIVE, N.O.S.

(<10% Isopropanol/Hydroxylamine Hydrochloride Solution)

**Hazard Class:** 3

**Subsidiary Risk:** 8

**ID Number:** UN2924

**Packing Group:** III

**I.M.O.:**

**Proper Shipping Name:** FLAMMABLE LIQUID, CORROSIVE, N.O.S.

(<10% Isopropanol/Hydroxylamine Hydrochloride Solution)

**Hazard Class:** 3

**Subsidiary Risk:** 8

**ID Number:** UN2924

**Packing Group:** III

**Additional Information:** There is a possibility that this product could be contained in a reagent set or kit composed of various compatible dangerous goods. If the item is NOT in a set or kit, the classification given above applies. If the item IS part of a set or kit, the classification would change to the following: UN3316 Chemical Kit, Class 9, PG II or III. If the item is not regulated, the Chemical Kit classification does not apply.

---

## 15. REGULATORY INFORMATION

### **U.S. Federal Regulations:**

**O.S.H.A.:** This product meets the criteria for a hazardous substance as defined in the Hazard Communication Standard. (29 CFR 1910.1200)

### **E.P.A.:**

**S.A.R.A. Title III Section 311/312 Categorization (40 CFR 370):** Immediate (Acute) Health Hazard Fire Hazard

**S.A.R.A. Title III Section 313 (40 CFR 372):** This product does NOT contain any chemical subject to the reporting requirements of Section 313 of Title III of SARA.

--

**302 (EHS) TPQ (40 CFR 355):** Not applicable

**304 CERCLA RQ (40 CFR 302.4):** Not applicable

**304 EHS RQ (40 CFR 355):** Not applicable

**Clean Water Act (40 CFR 116.4):** Not applicable

**RCRA:** Contains RCRA regulated substances. See Section 13, EPA Waste ID Number.

### **State Regulations:**

**California Prop. 65:** No Prop. 65 listed chemicals are present in this product.

**Identification of Prop. 65 Ingredient(s):** None

**California Perchlorate Rule CCR Title 22 Chap 33:** Not applicable

**Trade Secret Registry:** Not applicable

### **National Inventories:**

**U.S. Inventory Status:** All ingredients in this product are listed on the TSCA 8(b) Inventory (40 CFR 710).

**CAS Number:** Not applicable

**Canadian Inventory Status:** All ingredients of this product are DSL Listed.

**EEC Inventory Status:** All ingredients used to make this product are listed on EINECS / ELINCS.

**Australian Inventory (AICS) Status:** All ingredients are listed.

**New Zealand Inventory (NZIoC) Status:** All components either listed or exempt.

**Korean Inventory (KECI) Status:** Not listed - exempt. Quantity < 100 kg per annum.

**Japan (ENCs) Inventory Status:** All components either listed or exempt.

**China (PRC) Inventory (MEP) Status:** All components either listed or exempt.

---

## 16. OTHER INFORMATION

**References:** 29 CFR 1900 - 1910 (Code of Federal Regulations - Labor). Air Contaminants, Federal Register, Vol. 54, No. 12. Thursday, January 19, 1989. pp. 2332-2983. TLV's Threshold Limit Values and Biological Exposure Indices for 1992-1993. American Conference of Governmental Industrial Hygienists, 1992. Fire Protection Guide on Hazardous Materials, 10th Ed. Quincy, MA: National Fire Protection Association, 1991. In-house information. Technical Judgment.

**Complete Text of H phrases referred to in Section 3:** H226 Flammable liquid and vapour. H290 May be corrosive to metals. H351 Suspected of causing cancer. H317 May cause an allergic skin reaction.

**Revision Summary:** . Substantial revision to comply with EU Reg 1272/2008, Reg 1907/2006 and UN GHS (ST/SG/AC.10/36/Add.3).

### **Date of MSDS Preparation:**

**Day:** 14

**Month:** July

**Year:** 2014

**MSDS Prepared:** MSDS prepared by Product Compliance Department extension 3350

**CCOHS Evaluation Note:** It is offered under exemption from WHMIS labeling as specified in the Controlled Products Regulation (CPR) Section 17. It is offered under the interim policy that was established by Health Canada permitting use of GHS-formatted safety data sheets in Canada prior to revision of CPR to GHS. This product has been classified and labeled in accordance with the requirements of GHS (ST/SG/AC.10/36/Add.3). This SDS has been prepared in accordance with the requirements of GHS (ST/SG/AC.10/36/Add.3).

---

**Legend:**

NA - Not Applicable	w/w - weight/weight
ND - Not Determined	w/v - weight/volume
NV - Not Available	v/v - volume/volume

**USER RESPONSIBILITY:** Each user should read and understand this information and incorporate it in individual site safety programs in accordance with applicable hazard communication standards and regulations.

**THE INFORMATION CONTAINED HEREIN IS BASED ON DATA CONSIDERED TO BE ACCURATE. HOWEVER, NO WARRANTY IS EXPRESSED OR IMPLIED REGARDING THE ACCURACY OF THESE DATA OR THE RESULTS TO BE OBTAINED FROM THE USE THEREOF.**

**HACH COMPANY ©2014**

World Headquarters  
Hach Company  
P.O.Box 389  
Loveland, CO USA 80539  
(970) 669-3050

MSDS No: M00582

# MATERIAL SAFETY DATA SHEET

---

## 1. CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

**Product Name:** Titrant Solution Hardness 3 0.015 M EDTA  
**Catalog Number:** 42632

Hach Company  
P.O.Box 389  
Loveland, CO USA 80539  
(970) 669-3050

Emergency Telephone Numbers:  
(Medical and Transportation)  
(303) 623-5716 24 Hour Service  
(515)232-2533 8am - 4pm CST

**MSDS Number:** M00582  
**Chemical Name:** Not applicable  
**CAS Number:** Not applicable  
**Additional CAS No. (for hydrated forms):** Not applicable  
**Chemical Formula:** Not applicable  
**Chemical Family:** Not applicable  
**Intended Use:** Laboratory Use Hardness determination

---

## 2. HAZARDS IDENTIFICATION

This mixture is not classified as hazardous per GHS (UN publication ST/SG/AC.10/36/Add.3)

**GHS Classification:**

**Hazard categories:** .

**GHS Label Elements:**

Not applicable

**Hazard statements:** .

**Precautionary statements:** If skin irritation occurs: Get medical advice/attention.

**HMIS:**

**Health:** 0

**Flammability:** 0

**Reactivity:** 0

**Protective Equipment:** X - See protective equipment, Section 8.

**NFPA:**

**Health:** 0

**Flammability:** 0

**Reactivity:** 0

**Symbol:** Not applicable

**WHMIS Hazard Classification:** Not applicable

**WHMIS Symbols:** Not applicable

---

## 3. COMPOSITION / INFORMATION ON INGREDIENTS

**Hazardous Components according to GHS:**

**Ethylenediaminetetraacetic Acid, Disodium Salt**

**CAS Number:** 139-33-3

**Chemical Formula:**  $C_{10}H_{14}N_2Na_2O_8 \cdot 2H_2O$

**GHS Classification:** Acute Tox. 5-Orl, H303; Eye Irrit. 2A, H319; Aquatic Acute 2, H401

**Percent Range:** < 0.5

**Percent Range Units:** weight / weight



*PEL:* 15 mg/m<sup>3</sup> as total dust; 5 mg/m<sup>3</sup> as respirable dust  
*TLV:* 10 mg/m<sup>3</sup> as inhalable dust; 3 mg/m<sup>3</sup> as respirable dust

*WHMIS Symbols:* Not applicable

**EDTA Tetrasodium Salt**

*CAS Number:* 64-02-8  
*Chemical Formula:* C<sub>10</sub>H<sub>12</sub>N<sub>2</sub>Na<sub>4</sub>O<sub>8</sub> 2H<sub>2</sub>O  
*GHS Classification:* Acute Tox. 4-Orl, H302; Eye Dam. 1, H318  
*Percent Range:* < 0.5  
*Percent Range Units:* weight / weight  
*PEL:* 15 mg/m<sup>3</sup> as total dust; 5 mg/m<sup>3</sup> as respirable dust  
*TLV:* 10 mg/m<sup>3</sup> as inhalable dust; 3 mg/m<sup>3</sup> as respirable dust

*WHMIS Symbols:* Other Toxic Effects

*Hazardous Components according to GHS:* No  
**Demineralized Water**

*CAS Number:* 7732-18-5  
*Chemical Formula:* H<sub>2</sub>O  
*GHS Classification:* Non hazardous  
*Percent Range:* 65.0 - 75.0  
*Percent Range Units:* weight / weight  
*PEL:* Not established  
*TLV:* Not established

*WHMIS Symbols:* Not applicable

**Propylene Glycol**

*CAS Number:*  
*Chemical Formula:* C<sub>3</sub>H<sub>8</sub>O<sub>2</sub>  
*GHS Classification:* Not applicable  
*Percent Range:* 20.0 - 30.0  
*Percent Range Units:* weight / weight  
*PEL:* Not established  
*TLV:* Not established

*WHMIS Symbols:* Not applicable

---

#### 4. FIRST AID MEASURES

*General Information:* In the event of exposure, show this Material Safety Data Sheet and label (where possible) to a doctor.

*Advice to doctor:* Treat symptomatically.

*Eye Contact:* Flush eyes with water. Call physician if irritation develops.

*Skin Contact (First Aid):* Wash skin with plenty of water.

*Inhalation:* None required.

*Ingestion (First Aid):* Give large quantities of water. Call physician immediately.

---

#### 5. FIRE FIGHTING MEASURES

*Flammable Properties:* Material will not burn. Material will not burn.

*Fire Fighting Instruction:* As in any fire, wear self-contained breathing apparatus pressure-demand and full protective gear.

*Extinguishing Media:* Use media appropriate to surrounding fire conditions

*Extinguishing Media NOT To Be Used:* Not applicable

*Fire / Explosion Hazards:* This product will not burn or explode.

*Hazardous Combustion Products:* This material will not burn.

---

#### 6. ACCIDENTAL RELEASE MEASURES

**Spill Response Notice:**

Only persons properly qualified to respond to an emergency involving hazardous substances may respond to a spill according to federal regulations (OSHA 29 CFR 1910.120(a)(v)) and per your company's emergency response plan and guidelines/procedures. See Section 13, Special Instructions for disposal assistance. Outside of the US, only persons properly qualified according to state or local regulations should respond to a spill involving chemicals.

**Containment Technique:** Stop spilled material from being released to the environment.

**Clean-up Technique:** If permitted by regulation, Absorb spilled liquid with non-reactive sorbent material. Mark bag 'Non-hazardous trash', and dispose of as normal refuse. Decontaminate the area of the spill with a soap solution. Otherwise, Pick up spill for disposal and place in a closed container. Dispose of in accordance with local, state and federal regulations or laws.

**Evacuation Procedure:** Evacuate as needed to perform spill clean-up. If conditions warrant, increase the size of the evacuation.

**DOT Emergency Response Guide Number:** Not applicable

---

## 7. HANDLING AND STORAGE

**Handling:** Wash thoroughly after handling. Maintain general industrial hygiene practices when using this product.

**Storage:** Keep container tightly closed when not in use.

**Flammability Class:** Not applicable

---

## 8. EXPOSURE CONTROLS / PERSONAL PROTECTION

**Engineering Controls:** Maintain general industrial hygiene practices when using this product.

**Personal Protective Equipment:**

**Eye Protection:** safety glasses with top and side shields

**Skin Protection:** Not applicable

**Inhalation Protection:** adequate ventilation

**Precautionary Measures:** Wash thoroughly after handling.

**TLV:** Not established

**PEL:** Not established

**For Occupational Exposure Limits (OEL) for ingredients, see section 3 - Composition/Information on Ingredients.:**

---

## 9. PHYSICAL AND CHEMICAL PROPERTIES

**Appearance:** Clear, colorless liquid

**Physical State:** Liquid

**Molecular Weight:** Not applicable

**Odor:** None

**Odor Threshold:** Odorless

**pH:** 5.0

**Metal Corrosivity:**

**Corrosivity Classification:** Not classified as corrosive to metals according to GHS criteria.

**Steel:** Not determined

**Aluminum:** Not determined

**Specific Gravity/ Relative Density (water = 1; air = 1):** 1.026

**Viscosity:** Not determined

**Solubility:**

**Water:** Soluble

**Acid:** Soluble

**Other:** Not determined

**Partition Coefficient (n-octanol / water):** Not applicable

**Coefficient of Water / Oil:** Not applicable

**Melting Point:** Not determined

**Decomposition Temperature:** Not determined

**Boiling Point:** > 100° C (>212° F)

**Vapor Pressure:** Not determined

**Vapor Density (air = 1):** Not determined

**Evaporation Rate (water = 1):** 0.63

**Volatile Organic Compounds Content:** Not determined

**Flammable Properties:** Material will not burn.

**Flash Point:** Not applicable

**Method:** Not applicable  
**Flammability Limits:**  
**Lower Explosion Limits:** Not applicable  
**Upper Explosion Limits:** Not applicable  
**Autoignition Temperature:** Not applicable  
**Explosive Properties:**  
Not classified according to GHS criteria.  
**Oxidizing Properties:**  
Not classified according to GHS criteria.  
**Reactivity Properties:**  
Not classified as self-reactive, pyrophoric, self-heating or emitting flammable gases in contact with water according to GHS criteria.  
**Gas under Pressure:**  
Not classified according to GHS criteria.  
Not determined

---

## 10. STABILITY AND REACTIVITY

**Chemical Stability:** Stable when stored under proper conditions.  
**Mechanical Impact:** None reported  
**Static Discharge:** None reported.  
**Reactivity / Incompatibility:** None reported  
**Hazardous Decomposition:** No hazardous decomposition products known.  
**Conditions to Avoid:** Not applicable

---

## 11. TOXICOLOGICAL INFORMATION

**Toxicokinetics, Metabolism and Distribution:** No information available for mixture.  
**Toxicologically Synergistic Products:** None reported  
**Acute Toxicity:** Acute Toxicity Estimate (ATE) - Calculated from Ingredient Toxicity Data  
ATE Oral Rat LD50 = 165868 mg/kg  
**Specific Target Organ Toxicity - Single Exposure (STOT-SE):** Based on classification principles, the classification criteria are not met.  
**Specific Target Organ Toxicity - Repeat Exposure (STOT-RE):** Based on classification principles, the classification criteria are not met.  
**Skin Corrosion/Irritation:** Based on classification principles, the classification criteria are not met. Irritating to skin.  
**Eye Damage:** Based on classification principles, the classification criteria are not met.  
**Sensitization:** Based on classification principles, the classification criteria are not met.  
**CMR Effects/Properties (carcinogenic, mutagenic or toxic to reproduction):** No germ cell mutagenicity, carcinogenicity or reproductive toxicity data found.  
This product does NOT contain any IARC listed chemicals.  
This product does NOT contain any NTP listed chemicals.  
This product does NOT contain any OSHA listed carcinogens.  
**Symptoms/Effects:**  
**Ingestion:** No Effects Anticipated  
**Inhalation:** No effects anticipated  
**Skin Absorption:** No effects anticipated  
**Chronic Effects:** No effects anticipated  
**Medical Conditions Aggravated:** None reported

---

## 12. ECOLOGICAL INFORMATION

**Product Ecological Information:** --  
No ecological data available for this product. Mobility in soil: No data available  
**Ingredient Ecological Information:** Propylene Glycol: Pimephales promelas LC50 = 51400 mg/L; Daphnia magna 48 h LC50 34400 mg/L; Selenastrum capricornutum 96 h 19000 mg/L. Ethylenediaminetetraacetic Acid Disodium Salt: 96 Hr LC50 bluegill = 159 mg/L (Static).  
CEPA Statement: Water: Persistent. Not bioaccumulative. Not inherently toxic to aquatic organisms.  
CEPA Statement: EDTA Tetrasodium Salt, Ethylenediaminetetraacetic Acid Disodium Salt, Propylene Glycol: Not persistent. Not bioaccumulative. Not inherently toxic to aquatic organisms.

---

## 13. DISPOSAL CONSIDERATIONS

**EPA Waste ID Number:** Not applicable

**Special Instructions (Disposal):** Open cold water tap completely, slowly pour the material to the drain. Allow cold water to run for 5 minutes to completely flush the system. Otherwise, Dispose of material in an E.P.A. approved hazardous waste facility.

**Empty Containers:** Rinse three times with an appropriate solvent. Dispose of empty container as normal trash.

**NOTICE (Disposal):** These disposal guidelines are based on federal regulations and may be superseded by more stringent state or local requirements. Please consult your local environmental regulators for more information. In Europe: Chemical and analysis solutions must be disposed of in compliance with the respective national regulations. Product packaging must be disposed of in compliance with the country-specific regulations or must be passed to a packaging return system.

---

## 14. TRANSPORT INFORMATION

### **D.O.T.:**

**D.O.T. Proper Shipping Name:** Not Currently Regulated

--

**Hazard Class:** NA

**Subsidiary Risk:** NA

**ID Number:** NA

**Packing Group:** NA

### **T.D.G.:**

**Proper Shipping Name:** Not Currently Regulated

--

**Hazard Class:** NA

**Subsidiary Risk:** NA

**UN Number/PIN:** NA

**Packing Group:** NA

### **I.C.A.O.:**

**I.C.A.O. Proper Shipping Name:** Not Currently Regulated

--

**Hazard Class:** NA

**Subsidiary Risk:** NA

**ID Number:** NA

**Packing Group:** NA

### **I.M.O.:**

**Proper Shipping Name:** Not Currently Regulated

--

**Hazard Class:** NA

**Subsidiary Risk:** NA

**ID Number:** NA

**Packing Group:** NA

**Additional Information:** There is a possibility that this product could be contained in a reagent set or kit composed of various compatible dangerous goods. If the item is NOT in a set or kit, the classification given above applies. If the item IS part of a set or kit, the classification would change to the following: UN3316 Chemical Kit, Class 9, PG II or III. If the item is not regulated, the Chemical Kit classification does not apply.

---

## 15. REGULATORY INFORMATION

### **U.S. Federal Regulations:**

**O.S.H.A.:** This product does not meet the criteria for a hazardous substance as defined in the Hazard Communication Standard. (29 CFR 1910.1200)

### **E.P.A.:**

**S.A.R.A. Title III Section 311/312 Categorization (40 CFR 370):** This product is not hazardous under 29 CFR.1910.1200 and therefore is not covered by Title III under SARA.

**S.A.R.A. Title III Section 313 (40 CFR 372):** This product does NOT contain any chemical subject to the reporting requirements of Section 313 of Title III of SARA.

--

**302 (EHS) TPQ (40 CFR 355):** Not applicable

**304 CERCLA RQ (40 CFR 302.4):** Not applicable

**304 EHS RQ (40 CFR 355):** Not applicable

*Clean Water Act (40 CFR 116.4):* Not applicable

*RCRA:* Contains no RCRA regulated substances.

**State Regulations:**

*California Prop. 65:* No Prop. 65 listed chemicals are present in this product.

*Identification of Prop. 65 Ingredient(s):* None

*California Perchlorate Rule CCR Title 22 Chap 33:* Not applicable

*Trade Secret Registry:* Not applicable

**National Inventories:**

*U.S. Inventory Status:* TSCA Listed: Yes

*CAS Number:* Not applicable

*Canadian Inventory Status:* DSL Listed: Yes

*EEC Inventory Status:* EINECS Listed: Yes

*Australian Inventory (AICS) Status:* All ingredients are listed.

*New Zealand Inventory (NZIoC) Status:* All components either listed or exempt.

*Korean Inventory (KECI) Status:* All components of this product are either listed, listed as the anhydrous compound or exempt.

*Japan (ENCs) Inventory Status:* All components either listed or exempt.

*China (PRC) Inventory (MEP) Status:* All components either listed or exempt.

---

## 16. OTHER INFORMATION

*References:* Air Contaminants, Federal Register, Vol. 54, No. 12. Thursday, January 19, 1989. pp. 2332-2983. TLV's Threshold Limit Values and Biological Exposure Indices for 1992-1993. American Conference of Governmental Industrial Hygienists, 1992. In-house information. Vendor Information. Technical Judgment.

*Complete Text of H phrases referred to in Section 3:* H302 Harmful if swallowed. . H318 Causes serious eye damage. H319 Causes serious eye irritation.

*Revision Summary:* Substantial revision to comply with EU Reg 1272/2008, Reg 1907/2006 and UN GHS (ST/SG/AC.10/36/Add.3).

**Date of MSDS Preparation:**

*Day:* 06

*Month:* June

*Year:* 2014

*MSDS Prepared:* MSDS prepared by Product Compliance Department extension 3350

*CCOHS Evaluation Note:* It is offered under exemption from WHMIS labeling as specified in the Controlled Products Regulation (CPR) Section 17. It is offered under the interim policy that was established by Health Canada permitting use of GHS-formatted safety data sheets in Canada prior to revision of CPR to GHS. This product has been classified and labeled in accordance with the requirements of GHS (ST/SG/AC.10/36/Add.3). This SDS has been prepared in accordance with the requirements of GHS (ST/SG/AC.10/36/Add.3).

---

**Legend:**

NA - Not Applicable	w/w - weight/weight
ND - Not Determined	w/v - weight/volume
NV - Not Available	v/v - volume/volume

**USER RESPONSIBILITY:** Each user should read and understand this information and incorporate it in individual site safety programs in accordance with applicable hazard communication standards and regulations.

**THE INFORMATION CONTAINED HEREIN IS BASED ON DATA CONSIDERED TO BE ACCURATE. HOWEVER, NO WARRANTY IS EXPRESSED OR IMPLIED REGARDING THE ACCURACY OF THESE DATA OR THE RESULTS TO BE OBTAINED FROM THE USE THEREOF.**

HACH COMPANY ©2015