



Material	94898-540
Material description	2,2,4-Trimethylpentane
Grade	ACS
CAS Number	540-84-1
Molecular formula	(CH <sub>3</sub> ) <sub>2</sub> CHCH <sub>2</sub> C(CH <sub>3</sub> ) <sub>3</sub>
Molecular mass	114.23

## Characteristics

## Specifications

Assay (GC)	Min. 99.0 %
Color (Pt-Co)	Max. 10 APHA
Water-soluble titrable acid	Max. 0.0003 meq/g
Residue after evaporation	Max. 0.001 %
Sulphur compounds (as S)	Max. 0.005 %

## Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Michel Boisvert, Head of the Laboratory - Lachine  
VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3,  
QC, Canada

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.