

JACKETED BENCH SCALE REACTORS

Ace Glass Scale-Up Series™ reactors enable scaling up from 100mL to 6000mL on the bench and enjoy the same geometry when replicating those results in our Kilo Scale reactors. The reactor base system includes all of the essential items needed to begin building a customized reactor specifically for your application. For versatility, the system has been designed with a working temperature range of -60°C through 200°C, and wetted surfaces of borosilicate glass and PTFE.

Designed for diversity and ease of use. A simple base system building platform which allows any reactor system to be customized using catalog or custom designed parts.


Jacketed Base Systems include:

- Scale-Up Series™ Reactor Stand
- Jacketed Flask
- Rod Mount Clamps and Support Ring
- Quick Release Clamp and O-Ring
- Swing Latch Clamp
- Bottom Drain Valve

Base systems by desired capacity. (Rod Mount Motor only)

Capacity, mL	Stand Profile, D x W x H	Order Code
100	19.5" x 24.75" x 38"	10098-882
250	19.5" x 24.75" x 38"	10098-884
500	19.5" x 24.75" x 38"	10098-886
1000	19.5" x 24.75" x 38"	10098-888
2000	19.5" x 24.75" x 38"	10098-890
3000	19.5" x 24.75" x 48"	10098-892
4000	19.5" x 24.75" x 48"	10098-894
5000	19.5" x 24.75" x 48"	10098-896
6000	19.5" x 24.75" x 48"	30621-088

100mL to 6000mL



UNJACKETED BENCH SCALE REACTORS

Designed for diversity and ease of use. A simple base system building platform which allows any reactor system to be customized using catalog or custom designed parts.


Unjacketed Base Systems include:

- Scale-Up Series™ Reactor Stand
- Unjacketed Flask
- Rod Mount Clamps and Support Ring
- Swing Latch Clamp
- Bottom Drain Valve

Base systems by desired capacity. (Rod Mount Motor only)

Capacity, mL	Stand Profile, D x W x H	Order Code
100	19.5" x 24.75" x 38"	75877-608
250	19.5" x 24.75" x 38"	75877-618
500	19.5" x 24.75" x 38"	75877-620
1000	19.5" x 24.75" x 38"	75877-622
2000	19.5" x 24.75" x 38"	75877-624
3000	19.5" x 24.75" x 48"	75877-626
4000	19.5" x 24.75" x 48"	75877-628
5000	19.5" x 24.75" x 48"	75877-630
6000	19.5" x 24.75" x 48"	75877-632

100mL to 6000mL



Contact us for...

Custom glass and repair services.

avantor™
delivered by **vwr™**

VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada unless otherwise noted, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2022 Avantor, Inc. All rights reserved.