



Item Number P1314
Item Potassium Gluconate, Anhydrous, FCC
CAS Number 299-27-4
Molecular Formula $C_6H_{11}KO_7$
Molecular Weight 234.25
Assay ($C_6H_{11}KO_7$; Dried Basis).....Min. 98.0%

MAXIMUM LIMITS

Lead..... 2 mg/kg
Loss on Drying..... 3.0%
Reducing Substances..... 1.0%

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: **East Coast Facility:**
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843