



**Item Number** M1366  
**Item** Magnesium Sulfate, Heptahydrate, FCC  
**CAS Number** 10034-99-8  
**Molecular Formula**  $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$   
**Molecular Weight** 246.48

Assay ( $\text{MgSO}_4$ ; after Ignition).....Min. 99.5%  
Loss on Ignition..... 40.0-52.0%

MAXIMUM LIMITS

Lead..... 4 mg/kg  
Selenium..... 0.003%

***Spectrum Chemicals & Laboratory Products***

---

**Corporate Headquarters:**      **East Coast Facility:**  
14422 S. San Pedro St.      755 Jersey Ave.  
Gardena, CA 90248      New Brunswick, NJ 08901  
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843