



VWR ANALYTICAL

Product Specification

Catalog Number: BDH1103-4LG, BDH1103-4LP, BDH1103-19L, CABDH1103-4LG, CABDH1103-19L

Catalog Description: Acetonitrile

Grade: ACS

Date Current Specification Set: November 24 2009

All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.

Meets ACS Specifications for General Use

Characteristics	Specifications
Appearance:	Clear
Water (K.F.):	Less than 0.30%
Residue:	Less than 0.005% (0.0010g/50mL)
GC Purity:	Greater than 99.5%
Color (APHA):	Less than 10
Titration acid:	Less than 8 µeq/g (No more than 8.0 mL of 0.01N sodium hydroxide)
Titration base:	Less than 0.6 µeq/g (No more than 4.5 mL of 0.01N hydrochloric acid)

Signed on behalf of VWR:

Neal Fox
Quality Control Supervisor