

The most flexible tube rack designed for use with a multitude of tube diameters and heights

Holds virtually any brand tube 10-30 mm of varying heights

### Multiple Configurations

- Interchangeable tube rack inserts hold multiple tube sizes firmly in place to assist technicians constantly alternating between sample vessels
- Hinged shelves can be pivoted in either an upper or lower position to securely accommodate the size and height of tube being used

### Multi-Functional Design

- Open design allows for better filling visibility and tube graduations
- Rack wells have a unique teardrop design that enables better retention of tubes over traditional designs
- Rotate tube rack insert to an angle position to facilitate more comfortable pipetting and allow better visibility


**Autoclavable**



Lower position is ideal height when working with shorter microcentrifuge tubes



Ergonomic handles for ease of transporting the rack



Alphanumeric grid for quick sample identification and retrieval

Angled position for more comfortable pipetting

<p>Included: 1 Base + 6 Tube Rack Inserts</p>			
	Two 10-13 mm Tube Racks Fits 5 mL, 7.5 mL, 10 mL Tubes	Two 13-17 mm Tube Racks Fits 5 mL, 15 mL Tubes	Two 25-30 mm Tube Racks Fits 25 mL, 30 mL Tubes

## OneRack® Multi

US Patent No. 9776188

Item No.	Color	L x W x H in	L x W x H cm	# Inserts Included	Well Diam/Insert	Wells/Insert	UOM
4700330-002	Blue/Natural	12.2 x 4.2 x 3	30.9 x 10.7 x 7.7	6 (2 ea/ Diameter)	13 mm	20	1 ea
					17 mm	20	
					30 mm	6	
4700330-004	Purple/Natural	12.2 x 4.2 x 3	30.9 x 10.7 x 7.7	6 (2 ea/ Diameter)	13 mm	20	1 ea
					17 mm	20	
					30 mm	6	