



VWR® GENERAL-PURPOSE LABORATORY LABELING TAPE

As per the manufacturer, the below mentioned product meets the following criteria:

North American Catalog No:	89097-920
Lot Number:	711223
Date of Manufacture:	July 29, 2020
Description:	VWR TAPE RAINBOW 1/2X500IN CS24
Product Shelf life:	18 months
Country of Origin:	USA

QC Testing

The following data are nominal values based on PSTC, ASTM and other standard tests. The data should not be considered specifications.

Adhesion to Steel (oz/in of width) 50 (13.7 N/25mm) PSTC-101, Quick Stick (oz/in of width) To Kraft 8 (2.2 N/25mm), To Steel 23 (6.3 N/25mm)- Tensile Strength (lbs/in of width) 31 (136 N/25mm) PSTC- 131- Elongation (% at break) 3 PSTC-131

-Total Thickness (mils) 6.3 (0.16mm)

Product Specifications

Material: Flat back paper with synthetic rubber adhesive

BSE/TSE: The above product does not contain materials infected with or containing either BSE (bovine spongiform encephalopathy) or TSE (transmittable spongiform encephalopathy) a.k.a. mad cow disease.

Non-Pyrogenic: N/A

DNase & RNase Free: N/A

Cytotoxicity: N/A

Latex Free: The above listed does not contain raw natural rubber latex (NRL) or dried natural rubber (DNR).

Signed:



Jamie Ethier
VP Global Quality
VWR, Part of Avantor

Date: October 16, 2020