



Product Specification

Catalog Number: BDH9266

Catalog Description: BDH POTASS PHOSPHAT DBSC

Grade: ACS

Date Current Specification Set: NOV. 2012

All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.

Characteristics

Specifications

Chloride	<= 0.003 %
Heavy Metals (as Pb)	<= 0.0005 %
Insolubles	<= 0.01 %
Iron	<= 0.001 %
Loss on Drying	<= 1.0 %
Nitrogen Compounds	<= 0.001 %
pH (5%, Water) @25C	8.5 - 9.6
Purity	>= 98.0 %
Sodium	<= 0.05 %
Sulfate	<= 0.005 %

Signed on behalf of VWR:

Robin Crocker
Quality Control Manager