



Item Number S1340
Item Sodium Molybdate, Dihydrate, Crystal,
CAS Number 10102-40-6
Molecular Formula Na₂MoO₄•2H₂O
Molecular Weight 241.95

Assay (Na₂MoO₄•2H₂O)..... 99.5-103.0%
pH of a 5% Solution @ 25.0°C..... 7.0-10.5
MAXIMUM LIMITS
Insoluble Matter..... 0.005%
Chloride (Cl)..... 0.005%
Phosphate (PO₄)..... 5 ppm
Sulfate (SO₄)..... 0.015%
Ammonium (NH₄)..... 0.001%
Heavy Metals (as Pb)..... 5 ppm
Iron (Fe)..... 0.001%

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: **East Coast Facility:**
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843