

MSDS

PHOSPHATE STANDARDS, < 7000 ppm Phosphate (PO4), Orthophosphate Standards

Material Safety Data Sheet

Section 1: Chemical Product and Company Identification

Catalog Number: 5826, 5826.1, 5828, 5833, 5838, 5839, M0011204, P-206, P068100, P068150, R5826110, R5826300, R5827000, R5835000, R5835012, R5835025, R5835050, R5835125, R5835200, R5835250, R5835375, R5835500, SP068100, SP068120, SP068150	
Product Identity: PHOSPHATE STANDARDS, < 7000 ppm Phosphate (PO4), Orthophosphate Standards	
Manufacturer's Name: RICCA CHEMICAL COMPANY LLC	Emergency Contact(24 hr) -- CHEMTREC® Domestic: 800-424-9300 International: 703-527-3887
CAGE Code: 4TCW6, 0V553, 4XZQ2	
Address: 448 West Fork Dr Arlington, TX 76012	Telephone Number For Information: 817-461-5601
Date Prepared: 1/18/99	Revision: 13 Last Revised: 06/17/2008 Date Printed: 12/22/2009 12:49:44 pm

Section 2. Composition/Information on Ingredients

Component	CAS Registry #	Concentration	ACGIH TLV	OSHA PEL
Water, Deionized	7732-18-5	Balance	Not Available	Not Available
			Not Available	Not Available
Potassium Phosphate, Monobasic	7778-77-0	<1	Not Available	Not Available
			Not Available	Not Available

Section 3: Hazard Identification

Emergency Overview: Non-flammable, non-toxic, non-corrosive. Does not present any significant health hazards at these low concentrations. Wash areas of contact with water.

Target Organs: eyes, skin

Eye Contact: May cause slight irritation.

Inhalation: May cause allergic respiratory reaction to those allergic to phosphates.

Skin Contact: May cause slight irritation to those allergic to phosphates.

Ingestion: Ingestion of large quantity may cause stomach upset.

Chronic Effects/Carcinogenicity: None

IARC - No.

NTP - No.

OSHA - No.

Reproductive Information: Not Applicable.

Teratology (Birth Defect) Information: Not Applicable.

MSDS

PHOSPHATE STANDARDS, < 7000 ppm Phosphate (PO₄), Orthophosphate Standards

Section 4: First Aid Measures - In all cases, seek qualified evaluation.

Eye Contact: Irrigate immediately with large quantity of water for at least 15 minutes. Call a physician if irritation develops.

Inhalation: Remove to fresh air. Give artificial respiration if necessary. If breathing is difficult, give oxygen.

Skin Contact: Flush with plenty of water for at least 15 minutes. Call a physician if irritation develops.

Ingestion: Dilute with water or milk. Do not induce vomiting. Call a physician if necessary.

Section 5: Fire Fighting Measures

Flash Point: Not Available.

Method Used: Not Available.

LFL: Not Available.

UFL: Not Available.

Extinguishing Media: Use any means suitable for extinguishing surrounding fire.

Fire & Explosion Hazards: Not considered to be a fire or explosion hazard.

Fire Fighting Instructions: Use normal procedures/instructions.

Fire Fighting Equipment: Use protective clothing and breathing equipment appropriate for the surrounding fire.

Section 6: Accidental Release Measures

Absorb with suitable material (vermiculite, clay, etc.) and dispose of in accordance with local regulations.

Check with local agencies for the proper disposal of phosphate containing solutions.

Section 7: Handling and Storage

As with all chemicals, wash hands thoroughly after handling. Avoid contact with eyes and skin. Protect from freezing and physical damage.

Safety Storage Code: General

Section 8: Exposure Control/Personal Protection

Engineering Controls: No specific controls are needed. Normal room ventilation is adequate.

Respiratory Protection: Normal room ventilation is adequate.

Skin Protection: Chemical resistant gloves.

Eye Protection: Safety glasses or goggles.

Section 9: Physical and Chemical Properties

Appearance: Clear, colorless liquid

Odor: Odorless

Solubility in Water: Infinite

Specific Gravity: Approximately 1

pH: Not Available.

Boiling Point(°C): Approximately 100

Melting Point(°C): Approximately 0

Vapor Pressure: Not Applicable.

Section 10: Stability and Reactivity

Chemical Stability: Stable under normal conditions of use and storage.

Incompatibility: None identified.

Hazardous Decomposition Products: None identified.

Hazardous Polymerization: Will not occur.

Section 11. Toxicological Information

No LD50 (oral, rat) information reported in "Registry of Toxic Effects of Chemical Substances."

Section 12. Ecological Information

Ecotoxicological Information: No information found.

Chemical Fate Information: No information found.

MSDS

PHOSPHATE STANDARDS, < 7000 ppm Phosphate (PO₄), Orthophosphate Standards

Section 13. Disposal Considerations

Dilute with water and flush to sewer if regulations allow for the flushing of phosphate containing solutions. If not allowed, save for recovery or recycling in an approved waste disposal facility. Always dispose of in accordance with local, state and federal regulations.

Section 14. Transport Information

Part Numbers:

This product is not regulated.

Section 15. Regulatory Information (Not meant to be all inclusive - selected regulation represented)

OSHA Status:The above items either do not contain any specifically hazardous material or the potentially hazardous material is present in such low concentration that the items do not present any immediate threat to health and safety. These items do not meet the OSHA Hazard Communication Standard (29 CFR 1910.1200) definition of a hazardous material.

TSCA Status:All components of this solution are listed on the TSCA Inventory or are mixtures (hydrates) of items listed on the TSCA Inventory.

Sara Title III:

Section 302 Extremely Hazardous Substances:Not Applicable.

Section 311/312 Hazardous Categories:No

Section 313 Toxic Chemicals:Not Applicable.

California: None Reported.

Pennsylvania: None Reported.

RCRA Status: Not Applicable.

CERCLA Reportable Quantity: None Reported.

WHMIS: Not Applicable.

Section 16. Other Information

NFPA Ratings:

Health: 0 Flammability: 0 Reactivity: 0 Special Notice Key:None

HMIS Ratings:

Health: 0 Flammability: 0 Reactivity: 0 Protective Equipment:B (Protective Eyewear, Gloves)

Rev 1, 8-3-2000: Reformatted from WordPerfect® to Microsoft Word®; (Section 1) Revised emergency telephone number to CHEMTREC® 800-424-9300, revised concentration from 0 - 0.14 to <1%, changed title from 2 - 1000 ppm to < 7000 ppm; (Section 7) added storage code.

Rev 2, 10-09-2001: Reformatted to electronic data format.

Rev 3, 07-08-2002: (Section 1) added catalog number 5826.1.

Rev 4, 03-04-2003: (Section 1) added catalog number 5833.

Rev 5, 03-28-2003: (Section 1) added Solutions Plus catalog number P068150.

Rev 6, 11-13-2003: (Section 1) added Solutions Plus catalog number P068100.

Rev 7, 10-12-2004: (Section 1) added catalog number R5827000.

Rev 8, 02-11-2005: (Section 1) added SpectroPure catalog number SP068120.

Rev 9, 07-07-2005: (Section 1) added catalog numbers R5835000, R5835012, R5835025, R5835125, R5835200, R5835250, R5835375 and R5835500.

Rev 10, 09-07-2006: (Section 1) added Red Bird catalog number P-206.

Rev 11, 06-07-2007: (Section 1) added internal part number M0011204.

Rev 12, 06-14-2007: (Section 1) added catalog numbers R5826110 and R5826300.

Rev 13, 06-18-2008: (Section 1) added catalog number R5835050.



MSDS

For RICCA, SpectroPure, Red Bird, and Solutions Plus Brands
Emergency Contact(24 hr) -- CHEMTREC®
Domestic: 800-424-9300
International: 703-527-3887

PHOSPHATE STANDARDS, < 7000 ppm Phosphate (PO₄), Orthophosphate Standards

When handled properly by qualified personnel, the product described herein does not present a significant health or safety hazard. Alteration of its characteristics by concentration, evaporation, addition of other substances, or other means may present hazards not specifically addressed herein and which must be evaluated by the user. The information furnished herein is believed to be accurate and represents the best data currently available to us. No warranty, expressed or implied, is made and RICCA CHEMICAL COMPANY assumes no legal responsibility or liability whatsoever resulting from its use.