



Item Number CA185
Item Camphor, Synthetic, USP
CAS Number 76-22-2
Molecular Formula $C_{10}H_{16}O$
Molecular Weight 152.23

Melting Range..... 174°-179°C
Specific Rotation..... Optically inactive
Appearance of Solution..... To pass test

MAXIMUM LIMITS

Nonvolatile Residue..... 0.05%
Halogens..... 0.035%
Residual Solvents..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: **East Coast Facility:**
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843