

# Sterile 70% Isopropanol

Convenient trigger spray



## Applications

- Surface cleaning
- Use on gloved hands in sterile suites
- Wipedown for pass through to controlled environment
- Ideally suited for use with sterile cleanroom wipers



**North America**  
1210 South Park Drive  
Kernersville, NC 27284  
Tel 800.TEXTWIPE  
336.996.7046  
Fax 336.996.2297  
[www.texwipe.com](http://www.texwipe.com)  
[info@texwipe.com](mailto:info@texwipe.com)

**Europe/Middle East**  
Saffierlaan 5  
2132 VZ Hoofddorp  
The Netherlands  
Tel +31 88 1307 400  
Fax +31 88 1307 499

**Asia/Pacific**  
50 Tagore Lane  
#02-01 Entrepreneur Centre  
Singapore 787494  
Tel +65 6468 9433  
Fax +65 6468 6772

DS-3270 ©2011 ITW Texwipe Printed in USA  
Effective: September 2011

## Description

Sterile 70% Isopropanol contains 70% by volume USP-grade Isopropanol (isopropyl alcohol) and 30% USP-purified water. This solution is submicron filtered, filled into cleaned containers, double-bagged and gamma-irradiated. Each bottle is lot coded with an expiration date for easy recordkeeping.

## Non-atomizing adjustable spray

Customers can choose a stream delivery for precise application or the coarse spray to apply solution to larger areas. This is a fully assembled, ready-to-use applicator/container; no unpacking, handling or attachment of a separate spray head is required.

## Features

- Made with USP-grade Isopropanol
- Filtered through 0.2 micron filter
- Gamma-irradiated to a sterility level of  $10^{-6}$  with independent QC audits to ensure sterility
- Double-bagged in solvent-safe bags
- Non-atomizing adjustable spray
- Fully assembled, ready to use applicator/container
- Each bottle includes a set of 5 color indicator dots for easy identification
- Lot traceable
- Each lot tested for endotoxins
- Lot-specific information on each shipment simplifies recordkeeping:
  - » Certificate of Processing confirming radiation dosage
  - » Certificates of Compliance and Analysis
  - » Expiration date marked on each package

## Benefits

- Removes many surface contaminants
- Evaporates and leaves extremely low residues
- Prefiltered and sterilized for pass through and use in critical environments
- Eliminates mixing, filtration, sterilization, processing and QC expense of in-house production of sterile alcohol/water blends
- Ready to use sprayer eliminates handling and assembly of sprayhead and maintains sterility
- Meets FDA guideline for medical device bacterial endotoxin limit of not more than 20 EU/device

## Products

| Number        | Description  | Packaging           |
|---------------|--|---------------------|
| <b>TX3270</b> | Sterile 70% Isopropanol<br>16 fluid ounces (473 mL)<br>trigger spray | 12 polybottles/case |
| <b>TX3273</b> | Sterile 70% Isopropanol<br>32 fluid ounces (946 mL)<br>trigger spray | 12 polybottles/case |
| <b>TX8270</b> | Sterile 70% Isopropanol<br>8 fluid ounces (237 mL)<br>trigger spray  | 12 polybottles/case |

# Sterile 70% Isopropanol

TX3270, TX3273, TX8270

## Product Specifications

| Property         | Description                                     |
|------------------|---|
| IPA Quality      | Meets the requirements of USP Isopropyl Alcohol |
| Water Quality    | Meets the requirements of USP Purified Water    |
| Appearance       | Clear,colorless liquid                          |
| Odor             | Characteristic alcohol                          |
| % IPA Content    | 68%–72%   |
| Specific Gravity | 0.872–0.883                                     |
| Acidity          | Not more than 1.0mL of 0.020 N Sodium Hydroxide |
| Limit of NVR     | Does not exceed 5 mg (0.01%)                    |

**Certificate of Analysis provided with every lot.**

*ITW Texwipe is the only wiper company to be ISO 9001:2000,13485:2003,14001:2004 and OHSAS 18001:1999 registered.*