



Material	BDH1151-4LG
Material description	Toluene
Grade	ACS
CAS Number	108-88-3
Molecular formula	C <sub>6</sub> H <sub>5</sub> CH <sub>3</sub>
Molecular mass	92.14

Characteristics	Specifications
Assay (GC, corrected for water)	Min. 99.5 %
Substances darkened by sulphuric acid	Passes test
Color (Pt-Co)	Max. 10 APHA
Residue after evaporation	Max. 0.001 %
Sulphur compounds (as S)	Max. 0.003 %
Water	Max. 0.03 %

## Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Michel Boisvert, Head of the Laboratory - Lachine  
 VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3,  
 QC, Canada



For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.