



# Reagent Alcohol, Absolute AR (ACS)

Product No. 7019  
Specifications current as of: Jun 8  
2011

TEST	SPECIFICATION
Meets A.C.S. Specifications	
Assay:	
Specially Denatured Alcohol Formula 3A (200 proof)(v/v)	94.0 - 96.0 %
Methanol (CH <sub>3</sub> OH)	3.5 - 5.5 %
Isopropyl Alcohol (2-propanol)(v/v)	4.0 - 6.0 %
Color (APHA)	max. 10
Residue after Evaporation	max. 0.001 %
Specific Gravity at 15.56°C	0.782 - 0.793
Water (H <sub>2</sub> O)	max. 0.5 %

SDA Formula 3A consists of 5 volumes of methanol and 100 volumes of 200proof ethanol.

Product Information (not specifications):  
Appearance (clear, colorless liquid)

For Laboratory, Research or Manufacturing Use